

# Deleted Items for Prostheses List

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
--------------	--------------	---------------------	------	-------------------------

## AMO Australia Pty Ltd

### Ophthalmic

AO039	TECNIS\u00AE Toric 1-Piece IOL	1-piece IOL Posterior chamber Intraocular lens, incorporates a proprietary wavefront designed toric aspheric optic with a squared posterior optic edge designed to provide a 360 degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects.	+5.0 to +34.0 diopter powers in 0.5 diopter increments. Cylinder power: 1.00 diopter, 1.50 diopter, 2.25 diopter, 3.00 diopter, 3.75 diopter, 4.00 diopter, 4.50 diopter, 5.25 diopter and 6.00 diopter (as measured at the IOL plane). Overall diameter 13mm, Optic size 6mm.	Expanded into AO041, AO042
-------	--------------------------------	---	--	----------------------------

## Amplitude Australia

### Specialist Orthopaedic

LI077	Arcad osteosynthesis memory staples	Memory metal implants	8 to 10mm	Transferred to C1008
LI078	Arcad osteosynthesis memory staples	Memory metal implants	15 to 20mm	Transferred to C1006
LI080	Airlock Bone screw	Bone Screws Titanium Alloy (Locking screws)	Diameter 3-3.5mm, Length 8-30mm	Transferred to C1004
LI081	Airlock Bone screw	Bone Screws Titanium Alloy (non locking screws)	Diameter 3-3.5mm, Length 8-30mm	Transferred to C1007
LI082	Airlock Bone plate System	Osteosynthesis bone plate system	Left and right thickness 1.3mm to 1.35 mm, Length from 16mm to 55mm	Transferred to C1014
LI083	Lync intramedullary implant	Lync is a Titanium intramedullary implant	3.5mm to 3.9mm Length, angled 0 degrees and 10 degrees	Transferred to C1013
LI084	Nexis osteosynthesis screws	Bone screw titanium alloy (compressive)	Diameter 2.3mm, Length 10-30mm	Transferred to C1010
LI085	Nexis osteosynthesis screws	Bone screw titanium alloy (compressive)	Diameter 2.9-4 mm, Length 10-60mm	Transferred to C1009
LI086	Nexis osteosynthesis screws	Bone screw titanium alloy (headed)	Diameter 2.3mm, Length 10-30mm	Transferred to C1011
LI087	Nexis osteosynthesis screws	Bone screw titanium alloy (compressive)	Diameter 5-7 mm, Length 26-120mm	Transferred to C1012
LI088	Nexis osteosynthesis screws	Bone screw titanium alloy (headed)	Diameter 2.9-4mm, Length 10-60mm	Transferred to C1003
LI089	Nexis osteosynthesis screws	Bone screw titanium alloy (headed)	Diameter 5-7 mm, Length 26-120mm	Transferred to C1005
LI091	Arcad osteosynthesis memory staples	Stainless steel static staples	8 to 10 mm	Transferred to C1015

## Aurora BioScience Pty Ltd

### Neurosurgical

RA001	VNS Therapy Pulse/Pulse Duo	Model 102/102R Generator	One size only	Transferred to SA174
RA002	VNS Therapy Lead	Model 302	One size only	Transferred to SA175

## Austofix

### Hip

AX044	Hyperion Uncemented Stem	Titanium, gritblast, HA coated, collared, collarless and high offset	Sizes 9 to 20	Transferred to LI108
-------	--------------------------	--	---------------	----------------------

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
AX047	Hyperion Cobalt Chrome Femoral Head	Cobalt Chrome, 12/14 Taper	Sizes 22, 28 and 32 mm	Transferred to LH106
AX064	Austofix Ceramic System	Biolox ceramic femoral head	36mm	Transferred to LH104
AX065	Austofix Ceramic System	Biolox Ceramic Femoral Head	28 and 32 mm	Transferred to LH107
AX102	U2 Taper Sleeve	Taper Sleeve, Unipolar	-4, 0, +4, +8	Transferred to LH105

## Austspine Pty Ltd

### Specialist Orthopaedic

AJ005	Kasios TCP	TCP Bone graft substitute Granules	5cc	Transferred to LH561
AJ006	Kasios TCP	TCP Bone graft substitute Granules	10cc	Transferred to LH568
AJ007	JectOS	Injectable synthetic bone graft substitute	5cc	Transferred to LH560
AJ008	JectOS 10cc	Injectable synthetic bone graft substitute	10cc	Transferred to LH565
AJ010	KGBone	Synthetic bone substitute, custom-fit for KAGE cervical interbody cage	4mm-10mm	Transferred to LH559
AJ034	Jectos 20cc	Injectable synthetic bone graft substitute	20cc	Transferred to LH563
AJ041	i-FACTOR Putty Peptide Enhanced Bone Graft	P-15, ABM, Osteoinductive, Osteoconductive bone graft substitute	0 - 5.0mL	Transferred to LH562
AJ042	i-FACTOR Putty Peptide Enhanced Bone Graft	P-15, ABM, Osteoinductive, Osteoconductive bone graft substitute	5.0ml - 10ml	Transferred to LH566
AJ043	i-Factor Flex & Flex FR Peptide Enhanced Bone Graft \u2013 25mm	P-15, ABM, Osteoinductive, Osteoconductive bone graft substitute	25mm x 25mm x 4mm	Transferred to LH564
AJ044	i-Factor Flex & Flex FR Peptide Enhanced Bone Graft \u2013 25mm	P-15, ABM, Osteoinductive, Osteoconductive bone graft substitute	50mm x 25mm x 4mm, 100mm x 25mm x 4mm	Transferred to LH567

### Spinal

AJ009	Kasios KAGE	Cervical interbody cage, PEEK	Height 3mm - 10mm	Transferred to LH569
-------	-------------	-------------------------------	-------------------	----------------------

## B Braun Australia Pty Ltd

### Neurosurgical

BB109	Aneurysm Clip	Yasargil Phynox Mini	2.8mm x 110g to 7.0mm x 130g	Deleted
BB110	Aneurysm Clip	Yasargil Titanium	4.9mm x 150g to 20.0mm x 200g	Deleted
BB111	Aneurysm Clip	Yasargil Titanium Mini	2.8mm x 110g to 7.0mm x 110g	Deleted
BB112	Aneurysm Clip	Yasargil Reinforcing, Titanium/Phynox	One size only	Deleted
BB175	Yasargil Aneurysm Clip System	Phynox temporary clips are manufactured from phynox.	3mm - 20mm length, phynox temporary clips.	Deleted
BB176	Yasargil Aneurysm Clip System	Titanium temporary clips are manufactured from titanium	3mm- 25mm length, titanium temporary clips.	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
<b>Bausch &amp; Lomb (Australia) Pty Ltd</b>				
<b>Ophthalmic</b>				
BL025	SofPort Advanced Optic (AO) Violet Shield	Posterior chamber, aphakic, foldable silicone, anterior & posterior edge modification intraocular lens with aberration technology, a violet light	6mm optic 13mm overall length Modified C 5 degree angulation 0-30D, 0.5D steps	Deleted
BL033	FineVision POD F	Posterior chamber, UV absorbing, foldable, aspheric, acrylic, single piece lens with apodised diffractive, trifocal optic and blue light protection	An overall diameter of 11.40mm and optic diameter of 6.00mm. From 0 to +35 diopters with 0.5 diopters with 0.5 diopter increments.	Deleted
<b>Biotronik Australia Pty Ltd</b>				
<b>Cardiac</b>				
BT128	Lumax 740 HF-T	CRT-ICD with 40J output with Advanced Patient Management, Home Monitoring, automatic right and left ventricular Threshold Monitoring and conditionally MRI compatible	W x H x D (mm) 66 x 59 x 13, Volume: 39.8cc, Mass: 94g	Deleted
BT129	Lumax 740 VR-T DX	High Energy Single Chamber ICD with Atrial Diagnostics and Home Monitoring. Conditionally MRI compatible	W x H x D (mm) 66 x 55 x 13, Volume: 37.2cc, Mass: 92g	Deleted
BT130	Lumax 740 DR-T	Dual Chamber ICD with 40J output with Advanced Patient Management, Home Monitoring, Automatic Ventricular Threshold Monitoring and conditionally MRI compatible.	WxHxD (mm) 66 x 55 x 13, Volume: 37cc, Mass: 92g	Deleted
BT131	Lumax 740 VR-T	Single Chamber ICD with 40J output with Advanced Patient Management, Home Monitoring, Automatic Threshold Monitoring and MRI Conditional compatibility	WxHxD (mm) 66 x 55 x 13, Volume: 37cc, Mass: 92g	Deleted
<b>Boston Scientific Australia Pty Ltd</b>				
<b>Cardiac</b>				
BS150	Altrua 50 VDD	Dual chamber VDD pacemaker	10cc - 15cc	Deleted
BS288	SQ-RX Pulse Generator	Subcutaneous cardiac defibrillator	78.2mm x 65.5mm x 15.7mm	Deleted
<b>CareFusion Australia 316 Pty Ltd</b>				
<b>Specialist Orthopaedic</b>				
ZQ001	Ava-Tex Radiopaque Bone Cement Kit	13.8g Barium-Polymer Powder & 6ml Liquid Monomer (part No. BCPM003)	12mL	Transferred to HW579
ZQ002	AVAmox Plus	Radiopaque bone cement system	11G & 13G	Transferred to HW578
ZQ003	AVAmox	Radiopaque bone cement system	11G & 13G	Transferred to HW577

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
<b>Cochlear Limited</b>				
<b>Ear, Nose &amp; Throat</b>				
CO033	Baha Cordelle II sound processor	A powerful sound processor for the Baha hearing system with a bodyworn unit and a head worn transducer. Available in three colours. Body worn sound processor for the Baha hearing system.	One Size	Deleted
<b>Coloplast Pty Ltd</b>				
<b>Urogenital</b>				
CT014	NovaSilk Polypropylene Mesh	knitted monofilament polypropylene mesh	15cm x 15cm - One Size	Deleted
<b>ConMed Linvatec Australia Pty Ltd</b>				
<b>General Miscellaneous</b>				
LV127	Touchstone Linear Cutter Stapler 60mm	A device that incorporates a disposable staple cartridge that lays a double staggered row of staples on each side of a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	60mm length in staple line with choice of 2.5mm, 3.8mm and 4.5mm titanium staples	Deleted
LV128	Touchstone Linear Cutter Stapler 80mm	A device that incorporates a disposable staple cartridge that lays a double staggered row of staples on each side of a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	80mm length in staple line with choice of 3.8mm, 4.2mm and 4.5mm titanium staples	Deleted
LV129	Touchstone Linear Cutter Stapler 100mm	A device that incorporates a disposable staple cartridge that lays a double staggered row of staples on each side of a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	100mm length in staple line with choice of 3.8mm and 4.5mm titanium	Deleted
LV130	Touchstone Linear Cutter Stapler Cartridge (Reloads) 60mm	A disposable staple cartridge used in conjunction with a matched Touchstone Linear Cutter Stapler that lays a double staggered row of staples on each side a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	60mm length in staple line with choice of 2.5mm, 3.8mm and 4.5mm titanium staples	Deleted
LV131	Touchstone Linear Cutter Stapler Cartridge (Reloads) 80mm	A disposable staple cartridge used in conjunction with a matched Touchstone Linear Cutter Stapler that lays a double staggered row of staples on each side a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	80mm length in staple line with choice of 3.8mm, 4.2mm and 4.5mm titanium staples	Deleted
LV132	Touchstone Linear Cutter Stapler Cartridge (Reloads) 100mm	A disposable staple cartridge used in conjunction with a matched Touchstone Linear Cutter Stapler that lays a double staggered row of staples on each side a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	100mm length in staple line with choice of 3.8mm and 4.5mm titanium staples	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
LV133	Touchstone Linear Stapler Cartridge	A device housing a staple cartridge that places a double staggered row of staples into appropriate target tissue.	30, 45, 60 and 90mm staple lines in 2.5, 3.5 or 4.8mm staples	Deleted
LV134	Touchstone Linear Stapler	A device housing a staple cartridge that places a double staggered row of staples into appropriate target tissue.	30, 45, 60 and 90mm staple lines in 2.5, 3.5 or 4.8mm staples	Deleted
LV137	TOUCHSTONE Curved Linear Cutter	The Curved Linear Cutter is supplied sterile with a single-patient use cartridge with 4 staggered rows of titanium staples. The sizes are 20mm, 30mm, 37mm, 40mm, 43mm, 46mm to fit different tissues and surgical procedures. Among them, the CLC20, CLC25 and CLC30 are packaged with 2 inter-changeable cartridges	20mm-46mm	Deleted
LV138	Touchstone Curved Linear Cutter Reload	Reload for Touchstone Curved Linear Cutter	20mm-46mm	Deleted
LV140	Touchstone Endoscopic Linear Cutter	The Endoscopic Linear Cutter is supplied sterile with a single-patient use cartridge with 6 staggered rows of titanium staples. The sizes are 60mm, 160mm and 260mm which are designed for different surgical procedures.	60mm, 160mm, 260mm	Deleted
LV141	Touchstone Endoscopic Linear Cutter Reload	The Endoscopic Linear Cutter is supplied sterile with a single-patient use cartridge with 6 staggered rows of titanium staples. The sizes are 60mm, 160mm and 260mm which are designed for different surgical procedures.	45mm, 60mm	Deleted

## Control Medics Pty Ltd

### General Miscellaneous

CV005	Percutaneous Drainage Catheter Kit	Percutaneous Drainage Catheter Kit - Contains: Catheter with straightener, Metal stiffening Cannula, Dilator, Introducer Needle, Guide Wire, Flexible stiffening Cannula	6 French - 16 French and 19cm - 39cm Length	Deleted
-------	------------------------------------	--	---	---------

## CoreMed Pty Ptd

### Ear, Nose & Throat

RC005	Ventilation Tubes - Collar Button	Titanium w wire	1.00mm, 1.25mm, 1.50mm lumen	Deleted
RC008	Ventilation Tubes - Shepherd	Titanium with wire	1.00mm, 1.25mm lumen	Deleted

## Covidien Pty Ltd

### General Miscellaneous

AS128	Parietex-Multifilament Polyester Mesh	Parietex-Multifilament Polyester Anatomical Mesh Standard / With Slit and Folding / With Slit and No Touch System	13 x 09 cm, 15 x 10 cm, 16 x 12cm	Compressed into AS243
AS139	Parietex Composite 9FX	Dual sided polyester meshes with a resorbable adhesion-prevention hydrophilic film	9cm round	Compressed into AS242

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
AS140	Parietex Composite 12FX	Dual sided polyester meshes with a resorbable adhesion-prevention hydrophilic film	12cm round	Compressed into AS242
AS141	Parietex Composite 15FX	Dual sided Polyester meshes with a resorbable adhesion-prevention hydrophilic film	15cm round	Compressed into AS242
AS142	Parietex Composite 20FX	Dual sided polyester meshes with a resorbable adhesion-prevention hydrophilic film	20cm round	Compressed into AS245
AS144	Parietex Composite 2015FX	Dual sided polyester meshes with a resorbable adhesion-prevention hydrophilic film	20x15cm square	Compressed into AS245
AS145	Parietex Composite 2520FX	Dual sided polyester meshes with a resorbable adhesion-prevention hydrophilic film	25x20cm Oval	Compressed into AS244
AS146	Parietex Composite 3020FX	Dual sided polyester meshes with a resorbable adhesion-prevention hydrophilic film	30x20cm oval	Compressed into AS244
AS148	Parietex Composite 1510OSX	Dual sided polyester meshes with a resorbable adhesion-prevention hydrophilic film	15x10cm square	Compressed into AS242
AS149	Parietex Composite 2015OSX	Dual sided polyester meshes with a resorbable adhesion-prevention hydrophilic film	20x15cm square	Compressed into AS245
AS150	Parietex Composite 2520OSX	Dual sided polyester meshes with a resorbable adhesion-prevention hydrophilic film	25x20cm Square	Compressed into AS244
AS158	ABSORBATAACK \u2122	Laroscopic Absorbable Fixation Device	5mm, 15 tacks	Compressed into AS246
AS160	ABSORBATAACK \u2122	Laparoscopic Absorbable Fixation Device	5mm, 30 tacks	Compressed into AS246
AS166	Parietex Multifilament PolyesterMesh	Parietex-Multifilament Polyester Preshaped Mesh Three-dimensional/Ease Grip, with Slit	12 x 08 cms	Compressed into AS243
AS167	Parietex Multifilament Polyester Mesh	Parietex-Multifilament Polyester Mesh Three-dimensional/Ease Grip, with Slit	13 x 09 cms, 14 x 09 cms	Compressed into AS243
AS168	Parietex Multifilament Polyester Mesh	Parietex-Multifilament Polyester Folding Mesh/Bidimensional with integrated folding system with slit	14x10cms	Compressed into AS243
AS169	Parietex Multifilament Polyester Mesh Anatomical	Parietex-Multifilament Polyester Anatomical Mesh/Bidimensional/Three dimensional with integrated folding system.	14x10cms	Compressed into AS243

## Emergo Asia Pacific Pty Ltd

### Cardiothoracic

ER195	HeartMate II Left Ventricular Assist System - Implant Kit - Unsealed Grafts	Implant Kit with Unsealed Grafts	One size only	Transferred to SJ383
ER196	HeartMate II Left Ventricular Assist System - Implant Kit - Sealed Grafts	Implant Kit with Sealed Grafts	One size only	Transferred to SJ386
ER267	HeartMate II Left Ventricular Assist System - System Controller	System controller	One size only	Transferred to SJ384
ER268	HeartMate II Left Ventricular Assist System - Battery Set	Rechargeable HeartMate II-specific batteries	One size only, set of 4	Transferred to SJ385
ER269	HeartMate II Left Ventricular Assist System - Stabilization Belt	Fabric belt for Stabilizing percutaneous lead	Small / medium; set of 2 large; set of 2	Transferred to SJ381

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
ER270	HeartMate II Left Ventricular Assist System - Modular Belt	Modular belt to hold System Controller and Batteries	Small / Medium Large	Transferred to SJ382
ER271	HeartMate II Left Ventricular Assist System - Shower Bag	Shower Bag	One size only; set of 2	Transferred to SJ380

## Spinal

ER230	Spineway Mont Blanc Pedicle Screw System- Pedicle Screw, Polyaxial	Spineway Mont Blanc Pedicle Screw System- Pedicle Screw, Polyaxial	dia (4-8 mm) x length (25-110 mm)	Expanded into ER292, ER293
-------	--	--	-----------------------------------	----------------------------

## Gambro Pty Ltd

### Vascular

GA001	GamCath Dolphin Protect HighFlow Haemodialysis Catheter Kit	Haemodialysis catheter with a modified anti-thrombogenic surface (polymer system of polydimethylsiloxane blocks) and a non-release anti-microbial coating (Bismuth copolymer film) on internal and external surfaces of catheter.	109585 - 13Fr x 15cm double lumen; 109584 - 13FRx 20cm double lumen, 109705 - 13Fr x 25cm double lumen; 109706 - 13Fr x 15cm curved double lumen; 109707 - 13Fr x 20cm curved double lumen.	Transferred to BX339
GA002	GamCath High Flow Catheter Kit	GamCath High Flow Catheter is a polyurethane, kidney elliptical multi lumen, staggered tipped catheter kit, which ensures flows up to 430ml/min	11.5fr and 13fr guage. Available in 150mm, 175mm, 200mm and 250mm length and straight and curved configuration.	Transferred to BX338

## Global Manufacturing Technology Pty Ltd

### Hip

GM052	SMR Modular Femoral Stem Proximal Body	Titanium, grit blast, non HA,	40 to 70mm	Deleted
GM053	SMR Locking Screw	Locking Screw for Proximal Body to Stem	40 to 70mm	Deleted
GM055	SMR Modular Hip System	Calcar Component	40 to 70mm	Deleted
GM059	SMR Distal Locking Screw	Titanium, '4.5mm Cortical Bone Screws	20 to 80mm	Deleted
GM120	K2 Global Femoral Stem	Titanium, plasma spray, proximal HA	Size 0to size 12	Deleted
GM133	K2 Global Modular Femoral Stem	K2 Mid Stem	Size 2 to 12	Deleted
GM134	K2 Global Modular Femoral Stem	K2 Stem	K2 Stem Size 2, K2 Stem Size3, K2 Stem Size 4, K2 Stem Size 5, K2 Stem Size 6 , K2 Stem Size 7, K2 Stem Size 8, K2 Stem Size 10, K2 Stem Size 11,	Deleted
GM143	SMR Femoral Stem	SMR Distal Component 240mm	SMR distal component 14x 240mm, SMR distal component 16 x 240mm, SMR distal component 17x240mm, SMR distal component 18 x 240mm, SMR distal component 19 x 240mm, SMR distal component 20 x 240mm, SMR distal component 22 x 240mm	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
GM144	SMR Femoral Stem	SMR Distal Component <200mm	SMR Distal Component 14 x 160mm, SMR Distal Component 16 x 160mm, SMR Distal Component 17 x 160mm, SMR Distal Component 18 x 160mm, SMR Distal Component 19 x 160mm, SMR Distal Component 20 x 160mm, SMR Distal Component 22 x 160mm, SMR Distal Component 14 x 200mm, SMR Distal Component 16 x 200mm SMR Distal Component 18 x 200mm, SMR Distal Component 20 x 200mm, SMR Distal Component 22 x 200mm	Deleted
GM149	Absolut Cemented Femoral Stem	CoCr, highly polished, tapered, with neck	0 to 14	Deleted

### Specialist Orthopaedic

GM148	BonOs Inject	Injectable acrylic vertebroplasty cement	40gm	Deleted
-------	--------------	--	------	---------

### Spinal

GM153	SOCORE Spinal Spinal - Hooks, Pedicle/Laminar	Hooks, Pedicle/Laminar	Standard, Medium and Offset	Deleted
GM154	SOCORE Spinal System - Pedicle Screws	Pedicle Screw- Polyaxial. Pedicle Screw - Reduction	Diameter 4mm 25mm to 40mm Length. Diameter 5mm 30mm to 50mm Length. Diameter 6mm 30mm to 55mm length. Diameter 7mm 30mm to 50mm Length. Diameter 8mm 40mm	Deleted
GM155	SOCORE Spinal System - Crosslink	Crosslink	45mm, 55mm, 65mm	Deleted
GM158	SOCORE Spinal System- Locking Screw	Locking Screw	Standard & Mini	Deleted
GM159	SOCORE Spinal System - Monoaxial	Pedicle Screw - Monoaxial	Diameter 45mm 25mm to 40mm Length. Diameter 5mm 30mm to 50mm Length. Diameter 6mm 30mm to 55mm Length. Diameter 7mm 30mm to 50mm Length. Diameter 8mm 40mm	Deleted
GM160	SOCORE Spinal System - Standard and Pre-Bent Rod	Standard and Pre-Bent Rod	Diameter 5.5mm 40mm to 500mm Length.	Deleted

### Global Orthopaedic Technology

#### Hip

GO167	R120 Distal Centraliser	PMMA,	Various	Deleted
GO176	APEX K1tm Hip System	K1 Femoral Stem - titanium plasma uncemented non-HA	Stem std and lat sizes 1 to 12	Deleted

#### Knee

GO250	SHiVA Modular Knee Hinge Prosthesis - Femoral Component	Hinge Knee Femoral Component	Sizes 2,3 and 4 Left and right	Transferred to SF134
GO251	SHiVA Modular Knee Prosthesis - Tibial Tray	Tibial Tray	Size 2,3 and 4	Transferred to SF136
GO252	SHiVA Modular Knee Prosthesis - Rotating Bearing	Rotating Bearing	Sizes 2,3 and 4 x 10mm	Transferred to SF135
GO253	SHiVA Modular Knee Hinge Prosthesis - Patella Button	Patella Button	Sizes 2,3 and 4	Transferred to SF138



Billing Code	Product Name	Further Description	Size	Explanation of Deletion
GO254	SHiVA Modular Knee Hinge Prosthesis	SHiVA Modular Knee Hinge Prosthesis - Tibial Spacer	Size 2,3 and 4 x 5mm and 10mm thick	Transferred to SF132
GO255	Apex PS Revision System - Femoral Component	The femoral component of a total knee revisionsystem that incorprates a proportionally sized box	Sizes 1-6, plus Sizes 1+ to 5+ Left and Right	Expanded into GO286, GO287
GO263	SHiVA Hinge Knee Revision System - Stems	Knee revision Extension Stem	Straight 10 - 18mm Diameter; 50 to 150mm Length Offset 10 - 16 mm Diameter; 75mm Length	Transferred to SF137
GO264	SHiVA Hinge Knee Revision System - Spacer	SHiVA Hinge Knee Revision Spacer including locking screws	Size 2 -4, Thickness 5 and 10 mm for Distal Femoral and Posterior Femoral Reference Size 2-5 , Thickness 5 and 10 mm for Tibial Reference	Transferred to SF133

### Specialist Orthopaedic

GO092	C-Fix Cable Fixation System	Titanium Leader Cable Assembly	1-1.3mm	Deleted
-------	-----------------------------	--------------------------------	---------	---------

### Spinal

GO201	Cambria Interbody Cage	PEEK Cervical Interbody Cage	13mm x 12mm, 15mm x 13mm, 17mm x 13mm height - 5mm - 12mm	Deleted
GO202	Malibu Spinal System Polyaxial Screw	Polyaxial screw	4.5mm x 20mm - 8.5mm x 110mm in 5mm increments	Deleted
GO203	Malibu Spinal System Monoaxial Screw	Monoaxial Bone Screw	4.5mm x 25mm - 9.5mm x 80mm	Deleted
GO204	Malibu Spinal System Locking Cap	Locking cap with integral set screw	5.5mm and 6.35mm	Deleted
GO205	Malibu Spinal System Lateral Connector	Lateral Connector	30mm - 60mm	Deleted
GO206	Malibu Spinal System In-line Connector	In-line Connector	5.5mm and 6.35mm	Deleted
GO207	Malibu Spinal System Hooks	Hook	5.5mm x 5.5mm - 9.0mm x 9.5mm	Deleted
GO208	Malibu Spinal System Rods	Rod	5.5mm x 30mm - 600mm 6.35mm x 100mm	Deleted
GO209	Malibu Spinal System ESS CoNiCr Rods	Rod CoNiCr Precontored	5.5mm x 100mm - 600mm	Deleted
GO210	Malibu Spinal System CrossBar	CrossBar	5.5mm and 6.35mm diameter rods 23mm - 75mm length options	Deleted
GO211	Malibu Spinal System Polyaxial Reduction Screw	Malibu Spinal System Polyaxial or Uni Planar Reduction Screw	4.5mm x 25mm - 40mm in 5mm increments 5.5mm x 25mm - 50mm in 5mm increments 6.5mm x 35mm - 55mm in 5mm increments 7.5mm x 35mm - 55mm in 5mm increments	Deleted
GO212	Sierra Spinal System Occipital Plate	Occipital Plate	25mm-35mm 30mm-40mm 40mm-50mm 45mm-55mm 32mm-40mm	Deleted
GO213	Sierra Spinal System Set Screw	Set Screw	3.5mm	Deleted
GO214	Sierra Spinal System, Polyaxial screw	Polyaxial bone screw	3.5mm - 4.5mm diameter 10mm 50mm Length	Deleted
GO215	Sierra Spinal System Occipital Rod	Occipital Rod	60, 75, 90 degrees	Deleted
GO216	Sierra Spinal System Parallel Connector	In-line connector	3.5-3.5 - 6.35mm	Deleted
GO217	Sierra Spinal System Locking Cap	Locking Cap	3.5mm	Deleted
GO218	Sierra Spinal System Lateral Connector	lateral connector	3.5mm	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
GO219	Sierra Spinal System Rod	Precurved Rod	3.5mm x 20mm - 100m	Deleted
GO220	Sierra Spinal System Transitional Rod	Transitional Rod	3.5mm diameter / 5.5mm or 6.35mm	Deleted
GO221	Sierra Spinal System, Hooks	Hooks	4.0mm and 6.0mm	Deleted
GO222	Sierra Spinal System CrosBar	CrossBar	Small, medium, large, x-large	Deleted
GO223	Sierra Spinal System Occipital Screw	Occipital Screw	6mm x 14mm	Deleted
GO224	Sonoma Vaiaable Bone Screw	Vaiaable Bone Screw	4.0mm x 10mm - 20mm4.5mm x 10mm - 18mm	Deleted
GO225	Sonoma Limited Angle Bone Screw	Limited Angle Bone Screw	4.0mm x 10mm - 18mm4.5mm x 10mm - 18mm	Deleted
GO226	Sonoma Anterior Cervic Plate	Anterior Cervic Plate	<55	Deleted
GO227	Sonoma Anterior Cervical Plate	Anterior Cervical Plate	>55	Deleted
GO228	Pacifica Interbody Cage	PEEK Lumbar Interbody cages std and lordotic	Medium Lateral 9mm - 11mm x A/P 20mm - 24mm. Height 7mm - 17mm	Deleted
GO229	Ventura Interbody Cage	PEEK Lumbar Interbody Cage	Medial Lateral 9mm - 1mm x A/P 28mm. Height 7mm - 17mm	Deleted
GO230	Zuma bone screw	Bone screw	5.5mm - 6.0mm x 20mm - 30mm	Deleted
GO231	Zuma implant with integral plate	PEEK Interbody cage with integral plate	26mm x 24mm, 32mm x 25mm Height 10mm - 18mm	Deleted
GO232	Zuma Spinal System Locking Cover	Locking Cover	10mm - 20mm	Deleted
GO235	Sierra Spinal System Occipital Clip	Occipital Clip	One size only	Deleted
GO236	Redondo Spinal System	An ALIF Spacer manufactured from PEEK Optima	Width 26-55mm, Depth 18-24mm, Height 8-18mm 0-12 degrees lordosis	Deleted
GO237	Hollywood Spinal Posterior Interbody	A TLIF implant manufactured from PEEK Optima	11mm x 27mm x 7 to 18mm Lordotic in 1mm increments 11mm x 27mm x 7 70 18mm in 1mm increments	Deleted
GO238	Newport MIS Spinal System - Cap and Rod Assembly	Bulleted rod with an attached locking cap	Diameter 5.5mm, Length 30mm-110mm	Deleted
GO239	Newport MIS Spinal System - Cap	Locking Cap	One size only	Deleted
GO240	Newport MIS Spinal System - Double Lead	Double Lead Screw	Diameter 5.5-7.5mm Length 25-110mm	Deleted
GO241	Newport MIS Spinal System-MIS Screw	Cannulted, polyaxial, extended tab	Diameter 5.5-7.5mm Length 25-55mm	Deleted
GO242	Newport MIS Spinal System - Polyaxial Reduction Screws	Cannulated Polyaxial Reduction Screw	Diameter 5.5-7.5mm Length 25-110mm	Deleted
GO243	Newport MIS Spinal System - MIS Rod	Precontoured MIS Rod	30-130mm	Deleted
GO244	Newport MIS Spinal System - DLTA Polyaxial Screws	DLTA Cannulated Polyaxial Screw	Diameter 5.5-7.5mm, length 25-110mm	Deleted
GO248	Balboa Anterior Buttress Plate	An anterior intravertebral screw fixation device	Two sizes. 19mm x 15mm or 15mm x 24mm. Thickness tapers from 3.1mm at the top to 2.2mm at the bottom end of the plate.	Deleted
GO249	Balboa Anterior Buttress Plate Screws	Bone Screw	5.5mm x 20mm, 5.5mm x 25mm	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
GO265	Vu aPOD Prime Intervertebral Body Fusion Device- ALIF Cage	The Vu aPOD Prime is an intervertebral body (IBD) fusion device which can be used in combination with two titanium bone screws or a titanium SpinPlate	27x21mm, 31x24mm,35x27mm and 39x30mm Cage 8 degree, 12 degree and 15 degree lordis options 8mm,10mm,12mm,14mm,and 16mm height options	Deleted
GO266	Vu aPOD Prime Intervertebral Body Fusion Device- Locking Screw- Spinal	The Vu aPOD Prime locking screw is a bone screw used to fix the Vu aPOD Prime to the vertebral body doing away with the need for an anterior plate	5.5x20mm, 5.5x25mm,5.5x30mm,5.5x35mm,6.5x20mm,6.5x25mm,6.5x30mm,6.5x35mm Locking screw	Deleted
GO267	Zuma-C Implant	Integrated low-profile titanium anterior fixation plate with PEEK Optima Interbody spacer	13x12x16mm - 17x13x12mm	Deleted
GO268	Zuma-C Bone Screw	Variable angle bone screw	4.0-4.5mm x 10-18mm	Deleted
GO269	Zuma-C Locking Cover	Screw Locking Cover	6-12mm	Deleted
GO270	Vu ePOD Intervertebral Body Fusion Device	Peek Optima Lumbar Interbody Cage	8-16mm, non, 0, 5 and 10 degree lordosis	Deleted
GO271	Coral Spinal System - Fenestrated Pedicle Screw	Fenestrated Polyaxial Pedicle Screw	4.75mm - 9.5mm x 20-95mm	Deleted
GO272	Coral Spinal System - Locking Screw Assembly	Pedicle Locking Screw	One size fits all	Deleted
GO273	Coral Spine System Rod	Spinal Rod	5.5mm x 30-600mm	Deleted
GO274	Barricaid	Annular replacement prosthesis	Height: 14mm; Width 8mm, 10mm, 12mm	Transferred to ER305
GO278	Coral Spinal System-Cannulated Polyaxial Pedicle Screw	Cannulated Polyaxial Pedicle Screw	4.00mm - 9.5mm x 20-95mm	Deleted
GO279	Coral Spinal System-Pedicle & Lamina Hooks	Pedicle & Lamina Hooks	Small, Medium & Large	Deleted
GO280	Coral Spinal System-Solid Monoaxial Pedicle Screw	Solid Monoaxial Pedicle Screw	4.00mm - 9.5mm x 20-95mm	Deleted
GO281	Coral Spinal System-Solid Polyaxial Pedicle Screw	Solid Polyaxial Pedicle Screw	4.00mm - 9.5mm x 20-95mm	Deleted

## Guidant Australia Pty Ltd

### Cardiac

GU028	Guidant Aust Selute Picotip VDD Lead	VDD bipolar, passive fixation	.	Deleted
-------	--------------------------------------	-------------------------------	---	---------

## iNova Pharmaceuticals Australia Pty Ltd

### Ophthalmic

IL001	Crystalens HD	Accommodating lens with > 3 Dioptres of accommodation with dual mechanism of action: (1) forward axial movement of lens and (2) accommodative arching of optic. Posterior chamber, polyspheric optics, 3-piece, plate haptic, foldable UV light filtering, bic	10.0 to 33.0 Dioptres in .05 Dioptre steps	Deleted
IL002	Crystalens AO	Accommodating lens designed to provide > 3 Diopters of accommodation with dual mechanism of action: (1) forward axial movement of lens and (2) accommodative arching of optic. Posterior chamber, aspheric, aberration free optics, 3-piece, plate haptic, fol	10.0 to 33.0 D in 0.50 Diopter steps. 10.0 - 16.5 Diopters (12mm length) and 17.0 - 33.0 Diopters (11.5mm length).	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
<b>Insight Surgical Pty Ltd</b>				
<b>Ophthalmic</b>				
IS001	Preloaded Capsular Tension Ring	Capsular tension ring preloaded (acrylic)	13/11mm	Deleted
IS016	BiFlex 677ABY Hydrophilic Acrylic Foldable Aspheric Monofocal Lens with Blue light Filter	Posterior chamber IOL, 360° Special Square Edge, Aspheric Coma optimized	Optic diameter: 6.0 mm, Overall diameter: 13.0 mm Powers available: -10.0 - +9.0 D increment: 1.0 D +10.0 - +25.0 D increment: 0.5 D +25.5 - +30.0 D increment: 0.5 D +31.0-+45.0 D increment: 1.0 D	Deleted
IS017	Preloaded Monofocal Hydrophilic 677PY with Blue light Filter for Medjet PL-MA 2.2 single use injection system	Posterior chamber IOL, 360° Special Square Edge, Aspheric Coma optimized	Optic diameter: 6.0 mm, Overall diameter: 13.0 mm Powers available: -10.0 - +9.0 D increment: 1.0 D +10.0 - +25.0 D increment: 0.5 D +25.5 - +30.0 D increment: 0.5 D +31.0-+45.0 D increment: 1.0 D	Deleted
IS021	Medicontur QFlex 640PMY Hydrophilic Acrylic Foldable Multifocal Lens with Blue light Filter	Posterior chamber IOL, 360° Special Square Edge, Aspheric Coma optimized	Optic diameter: 6.0 mm, Overall diameter: 13.0 mm Powers available: -10.0 - +9.0 D increment: 1.0 D +10.0 - +25.0 D increment: 0.5 D +25.5 - +30.0 D increment: 0.5 D +31.0- +45.0 D increment: 1.0 D	Deleted
IS032	Medicontur Preloaded Monofocal Hydrophilic 640P without Blue light Filter for Medjet PL-MA 2.2 single use injection system	Posterior chamber IOL, 360° Special Square Edge, Aspheric Coma optimized	Optic diameter: 6.0 mm, Overall diameter: 13.0 mm Powers available: -10.0 - +9.0 D increment: 1.0 D +10.0 - +25.0 D increment: 0.5 D +25.5 - +30.0 D increment: 0.5 D +31.0-+45.0 D increment: 1.0 D	Deleted
IS033	BiFlex 877FABY Hydrophobic Acrylic Foldable Aspheric Monofocal Lens with Blue light Filter	Posterior Chamber IOL, special 360° sharp edge, Aspheric - Aberration Balanced	Optic diameter: 6.0 mm, Overall diameter: 13.0 mm Powers available: -0.0 - +9.0 D increment: 1.0 D +10.0 - +25.0 D increment: 0.5 D +25.5- +30.0 D increment: 0.5 D +31.0 - +35.0 D increment: 1	Deleted
IS035	BiFlex 877FAB Hydrophobic Acrylic Foldable Aspheric Monofocal Lens without Blue light Filter	Posterior Chamber IOL, special 360° sharp edge, Aspheric - Aberration Balanced	Optic diameter: 6.0 mm, Overall diameter: 13.0 mm Powers available: -0.0-+9.0 D increment: 1.0 D +10.0-+25.0 D increment: 0.5 D +25.5-+30.0 D increment: 0.5 D +31.0-+35.0 D increment: 1.0 D	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
<b>Integra Neurosciences Pty Ltd</b>				
<b>Specialist Orthopaedic</b>				
IG074	Tac-Pin (Threaded)	Compression pin-titanium	1.7mm diam 15-25mm length	Deleted
<b>Johnson &amp; Johnson Medical Pty Ltd</b>				
<b>General Miscellaneous</b>				
JJ891	Endopath Linear Cutter	Metal/plastic cutter-stapler with titanium staples	45mm; short shaft	Compressed into MN197
JJ893	Endopath Linear Cutter	Metal/plastic cutter-stapler with titanium staples	45mm; long shaft	Compressed into MN197
JJ897	Endopath Linear Cutter	Metal/plastic cutter-stapler with titanium staples	35mm	Compressed into MN197
MN091	PHYSIOMESH	Flexible composite mesh	7x15cm, rectangular shape	Deleted
MN092	PHYSIOMESH	Flexible Composite Mesh	10x15cm, oval shape	Deleted
MN093	PHYSIOMESH	Flexible Composite Mesh	15x15cm, square shape	Deleted
MN094	PHYSIOMESH	Flexible Composite Mesh	15x20cm, rectangular shape	Deleted
MN095	PHYSIOMESH	Flexible Composite Mesh	15x20cm, Oval Shape	Deleted
MN096	PHYSIOMESH	Flexible Composite Mesh	20x25cm, oval shape	Deleted
MN097	PHYSIOMESH	Flexible Composite Mesh	20x30cm, rectangular shape	Deleted
MN098	PHYSIOMESH	Flexible Composite Mesh	25x35cm, oval shape	Deleted
MN099	PHYSIOMESH	Flexible Composite Mesh	30x35cm, rectangular shape	Deleted
MN100	PHYSIOMESH	Flexible Composite Mesh	30x50cm, rectangular shape	Deleted
<b>Johnson &amp; Johnson Medical Pty Ltd t/a DePuy Synthes</b>				
<b>Specialist Orthopaedic</b>				
MN040	Bristow-Latarjet Instability Shoulder Set - Screw	Cannulated titanium locking screw	28mm to 48mm	Expanded into MN199, MN200
<b>Spinal</b>				
JJ835	Viper Spine System	Polyaxial and Uniplanar Screw, with and without tab	20-130mm	Compressed into MN198
JJ836	Viper Spine System	Polyaxial screw with cortical fixation	20-120mm	Compressed into MN198
MN156	Viper Spinal System	F2 Transfacet Pedicular Screw	5 to 7mm	Deleted
<b>Lifehealthcare Pty Ltd</b>				
<b>Cardiac</b>				
LH521	BioMatrix Flex Drug Eluting Coronary Stent System	Coronary Artery Stent, Drug Eluting	(8 - 38)mm length x (2.25 - 4.0)mm diameter	Transferred to HE001
LH522	Gazelle Coronary Stent	Stainless Steel	Diam 2.25-4.0mm Length 8-28mm	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
<b>Plastic and Reconstructive</b>				
LH493	Neuro Flap Fix Kit (Sterile)	Neuro Flapfix Kit (sterile) consists of three low-profile 2-hole plates & six self-drilling 1.5-mm maxDrive screws.	Plates: length 15mm / width 1mm /thick 0.6mm Screws: length 4mm/ width 1.5mm	Deleted
<b>Lima Orthopaedics Australia Pty Ltd</b>				
<b>Specialist Orthopaedic</b>				
LC049	DCD Cerclage Cable	The DCD cerclage cable is used to treat femoral fractures, humeral and tibial fractures, diaphyseal reconstruction in revision hip arthroplasty and diaphyseal reconstruction in revision shoulder arthroplasty. The spacers minimise contact with the periostium to reduce the potential of interruption to the periosteal blood supply to facilitate bone healing.	Single size only.	Deleted
LC056	Self Breaking Screw - Mini	Self Breaking Screw Mini	Diameter 2.7 mm x length 11-22 mm	Deleted
LC064	Hot Memory Staple 4 leg - Medium	Hot Memory Staple 4 leg - Medium	Width 11-20 mm	Deleted
LC065	ADS INTERPHALANX COLD MEMORY STAPLE	INTERPHALANX Cold Memory Staple	0-35 Degree	Deleted
LC067	DTS - Double threaded cannulated screw with guidewire	DTS - Double threaded cannulated screw with guidewire	Diameter 3, 3.5, x Length 12 to 40mm	Deleted
LC068	DTS - DOUBLE THREADED CANNULATED SCREW	DTS - DOUBLE THREADED - SELF TAPPING OR SELF DRILLING CANNULATED SCREW	Diameter 4.5 and 6.5mm x Length 30 to 100mm	Deleted
LC071	TNP\u00AE Cannulated Self drilling Screw with Ultra Flathead	TNP\u00AE Cannulated Self drilling Screw with Ultra Flathead	Diameter 4.5, 6mm x 20-120mm Length	Deleted
LC072	Ligament Fixation Staple - Small	Ligament Fixation Staple - Small	6 - 11 mm width	Deleted
LC073	Straight Staples - Small	Straight Staples - Small	7.5 to 10 mm width	Deleted
LC074	Osteotomy Offset Staple - Small	Osteotomy Offset Staple - Small	5 - 11 mm	Deleted
LC076	Steinmann Pin	Steinmann	Diameter 2 - 5mm x Length 150 - 250mm	Deleted
LC077	Self Breaking Screw - Micro	Self Breaking Screw - Micro	Dia 2.0 mm x Length 7-24 mm	Deleted
LC079	Mono-Cortical Hot Memory Staple - 2 leg large	Mono-Cortical Hot Memory Staple - 2 leg large	Width 20-30 mm	Deleted
LC080	Mono-cortical Hot Memory Staple - 2 leg Medium	Mono-cortical Hot Memory Staple - 2 leg Medium	Width 11-20 mm	Deleted
LC081	Mono-cortical Hot Memory Staple - 2 leg Small	Mono-cortical Hot Memory Staple - 2 leg Small	Width 8-11 mm	Deleted
LC082	Asymmetric Hot Memory Staple - Medium	Asymmetric Hot Memory Staple - 2 leg - Medium	Width 11-12mm	Deleted
LC083	Ligament Fixation Staple - Medium	Ligament Fixation Staple - Medium	11 - 16 mm width	Deleted
LC084	Straight Staples - Medium	Straight Staples - Medium	12.5 mm	Deleted
LC085	Osteotomy Offset Staple - Medium	Osteotomy Offset Staple - Medium	11 - 20 mm	Deleted
LC086	Osteotomy Offset Staple - Large	Osteotomy Offset Staple - Large	26 - 36 mm	Deleted
LC089	Kirschner Wire	Kirschner Wire - Various	Diameter 1 to 3mm x 80 to 310mm Length	Deleted
LC099	Osteosynthesis Bone Screws - Micro	Osteosynthesis Bone Screws - Micro	Dia 1.5/1.8/2.0mm x 4-40mm length	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
LC100	Osteosynthesis Bone Screws - Mini	Osteosynthesis Bone Screws - Mini	Dia 2.7mm x 6-40mm length	Deleted
LC101	Osteosynthesis Bone Screws Cancellous - Small	Osteosynthesis Bone Screws Cancellous - Mini	Dia 3.5mm x 10-120mm length	Deleted
LC102	Osteosynthesis Bone Screw Cortical - Small	Osteosynthesis Bone Screw Cortical - Small	Dia 3.5mm x 10-120mm length	Deleted
LC105	Supernail / Uninail Cephalic Screw	Supernail / Uninail Cephalic Screw	Diameter 12mm x Length 70 - 120mm	Deleted
LC106	Supernail / Uninail - Self Tapping Bone Screw	Supernail / Uninail - Self Tapping Bone Screw	Diameter 4.5 / 6mm x Length 24 - 100mm	Deleted
LC118	Superscrew - Mini	Cannulated Cortical Partial or Full Thread	Dia. 2.7mm x 10-30mm length	Deleted
LC119	Superscrew - small	Cannulated Cortical / Cancellous Lag Screw Partial or Full Thread	Dia. 3.5 / 4.0 mm x 10-70mm length	Deleted
LC120	Superscrew - Standard	Cannulated Cortical / Cancellous Lag Screw Partial or Full Thread	Dia. 4.5 / 5.0 / 6.5 / 7.2 mm x 25-140mm Length	Deleted
LC121	Versa Screw Washer	Versa Screw Washer	One Size	Deleted
LC122	Versa Screw	Self Drilling self tapping Cannulated Compression Screw	6.5 to 8mm Dia x 40 to 120mm length	Deleted
LC123	TwinLock Screw - Mini	TwinLock Screw - Mini	2.5/2.8mm Dia x 8 to 20 mm length	Deleted
LC126	R2 Subtalar Implant	R2 Subtalar Implant	Size 7/13mm to 12/18mm	Deleted

## LMT Surgical Pty Ltd

### Specialist Orthopaedic

LM069	Normed Foot screws	Titanium locking bone screws	2.7-4.0mm	Transferred to ZI941
LM127	Normed fracture fixation plates	Titanium locking foot plates	4/6 hole and 2.7mm screws	Transferred to ZI946
LM128	Normed fracture fixation plates	Titanium locking foot plates	4/6 hole and 3.5mm screws	Transferred to ZI945
LM129	Normed fracture fixation plates	Titanium locking foot plates	16 hole	Transferred to ZI942
LM130	Normed fracture fixation plates	Titanium locking foot plates	8/10/12 holes and 3.5mm screw	Transferred to ZI944
LM168	Ankle Fix - Plate \u2265 7 to \u2264 15 holes	Ankle Fusion plate	Standard and Large	Transferred to ZI948
LM169	Ankle Fix -Plate \u2265 16 holes	Ankle fusion plate	One size only	Transferred to ZI943
LM170	Screws - Standard	Standard Screws	10 - 70mm	Transferred to ZI940
LM171	Screws -Locking	Locking Screw	10 - 70mm	Transferred to ZI947
LM172	NEUMANN Wrist Arthrodesis System - Plate	Wrist arthrodesis plate	small and large	Deleted
LM173	Screws - Standard - Micro	2mm diameter standard screw	5 - 19mm	Deleted
LM174	Screws - Standard - Mini	2.7mm diameter standard screw	8 - 30mm	Deleted
LM175	Screws -Locking -mini	2.7mm diameter locking screw	8 - 30mm	Deleted
LM178	PIDI-KIT System - Plate	Proximal radial head and neck fracture plate	Small, medium and large	Deleted
LM179	Compression Screw - Small	3.9mm diameter cannulated compression screw	Length: 14 - 22mm	Deleted
LM180	Compression Screw - Mini	2.7mm diameter cannulated compression screw	Length:14 - 24mm	Deleted
LM272	OsteoMed FPS	Titanium locking footplates	3.5mm and 4mm	Expanded into LM304, LM305
LM273	OsteoMed FPS	Titanium Locking footplates	3.5mm and 4mm	Expanded into LM306, LM307, LM308, LM309

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
<b>Medacta Australia Pty Ltd</b>				
<b>Knee</b>				
MU044	Cin\u00E9tique Knee System Femoral Component	Femoral Component, cemented, CoCr, with TiNi (TiNibN) surface coating	Sizes 1 to 6, Left and Right	Deleted
<b>Med-Chem Surgical Pty Ltd</b>				
<b>General Miscellaneous</b>				
ML004	Touchstone Linear Cutter Stapler 60mm	A device that incorporates a disposable staple cartridge that lays a double staggered row of staples on each side of a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	60mm length in staple line with choice of 2.5mm, 3.8mm and 4.5mm titanium staples	Deleted
ML005	Touchstone Linear Cutter Stapler 80mm	A device that incorporates a disposable staple cartridge that lays a double staggered row of staples on each side of a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	80mm length in staple line with choice of 3.8mm and 4.5mm titanium staples	Deleted
ML006	Touchstone Linear Cutter Stapler 100mm	A device that incorporates a disposable staple cartridge that lays a double staggered row of staples on each side of a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	100mm length in staple line with choice of 3.8mm and 4.5mm titanium	Deleted
ML007	Touchstone Linear Cutter Stapler Cartridge (Reloads) 60mm	A disposable staple cartridge used in conjunction with a matched Touchstone Linear Cutter Stapler that lays a double staggered row of staples on each side a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	60mm length in staple line with choice of 2.5mm, 3.8mm and 4.5mm titanium staples	Deleted
ML008	Touchstone Linear Cutter Stapler Cartridge (Reloads) 80mm	A disposable staple cartridge used in conjunction with a matched Touchstone Linear Cutter Stapler that lays a double staggered row of staples on each side a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	80mm length in staple line with choice of 3.8mm and 4.5mm titanium staples	Deleted
ML009	Touchstone Linear Cutter Stapler Cartridge (Reloads) 100mm	A disposable staple cartridge used in conjunction with a matched Touchstone Linear Cutter Stapler that lays a double staggered row of staples on each side a cutting knife. Used for extensive bowel, gastric, gynaecological and thoracic surgical stapling procedures.	100mm length in staple line with choice of 3.8mm and 4.5mm titanium staples	Deleted
ML010	Touchstone Linear Stapler Cartridge	A device housing a staple cartridge that places a double staggered row of staples into appropriate target tissue.	30, 45, 60 and 90mm staple lines in 2.5, 3.5 or 4.8mm staples	Deleted
ML011	Touchstone Linear Stapler	A device housing a staple cartridge that places a double staggered row of staples into appropriate target tissue.	30, 45, 60 and 90mm staple lines in 2.5, 3.5 or 4.8mm staples	Deleted



Billing Code	Product Name	Further Description	Size	Explanation of Deletion
ML012	Touchstone Circular Stapler	A surgical implanting device for placing a specific size of double staggered row of titanium staples in a circular fashion in the large bowel, small bowel or oesophagus to maintain function of the digestive tract. A suitable size is used to match the lumen size of the patient's digestive tract	external dimension - 25mm, 29mm or 33mm	Deleted
ML013	Touchstone Prolapse Haemorrhoid Stapler	A surgical implanting device for placing a specific size of double staggered row of titanium staples in a circular fashion in the Anal and/or Rectal region for the circumferential excision of the prolapsed distal rectal mucosa, as described by Dr. Longo, Italy.	33mm is the external diameter of the device.	Deleted

## Medical Vision Australia Plastic and Cosmetic Pty Ltd

### Ophthalmic

MW002	Biopore Sphere	An implantable device made of high-density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device is available in sizes to meet the patient's specific needs. Used for orbital reconstruction following enucleation and evisceration procedure.	14mm, 16mm, 18mm, 20mm, 22mm, 24mm (diameter)	Transferred to S2005
-------	----------------	---	---	----------------------

### Plastic and Reconstructive

MW001	Biopore Flexible Block	An implantable device intended to support, contour, and provide structure for the craniofacial anatomy during reconstructive and plastic surgery procedures. Made of high density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device is available in various shapes and can be cut by the surgeon to meet the patient's specific needs.	38 x 50 x 3mm and 76 x 50 x 9mm	Transferred to S2003
MW003	Biopore Mandible (Left / Right)	An implantable device intended to support, contour, and provide structure for the mandible. Made of high-density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device can be cut by the surgeon to meet the patient's specific needs.	74 x 33 x 5mm	Transferred to S2009
MW004	Biopore Malar Implant	An implantable device intended to support, contour and provide structure for the face including the cheeks and eye sockets. Made of high-density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device is available in various shapes and can be cut by the surgeon to meet the patient's specific needs.	Left and Right (47mm x 24mm x 4mm)	Transferred to S2004

<b>Billing Code</b>	<b>Product Name</b>	<b>Further Description</b>	<b>Size</b>	<b>Explanation of Deletion</b>
MW005	Biopore Nasal Dorsum Shape	An implantable device intended to support, contour, and provide structure for the nose and zygoma. Made of high-density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device is available in various shapes and can be cut by the surgeon to meet the patient's specific needs.	53 x 5mm, 66 x 6 mm, 76 x 15mm, 40 x 9mm	Transferred to S2010
MW006	Biopore Ultra thin Sheet	An implantable device intended to support, contour, and provide structure for the craniofacial anatomy during reconstructive and plastic surgery procedures. Made of high-density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device is available in various shapes and can be cut by the surgeon to meet the patient's specific needs.	38 x 50mm	Transferred to S2002
MW007	Biopore Ultra Thin Sheet Large	An implantable device intended to support, contour, and provided structure for the craniofacial anatomy during reconstructive and plastic surgery procedures. Made of high- density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device is available in various shapes and can be cut by the surgeon to meet the patient's specific needs	76 x 127mm	Deleted
MW008	Biopore Orbital Rim, Extended	An implantable device intended to support, contour, and provide structure for the face including the cheeks and eye socket. Made of high-density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device is available in various shapes and can be cut by the surgeon to meet the patient's specific needs.	Left and Right (44mm x 38mm)	Transferred to S2006
MW009	Biopore Custom Made Implant	An implantable device made of high-density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device is custom made/biomodelled from individual patient's CT scan, on doctors prescription, used for correction and reconstruction of any facial, head or any non weight bearing bone structure.	Custom designed and biomodelled from patient's CT scan image.	Transferred to S2007
MW010	Biopore Chin Two Piece Anatomical contour.	An implantable device intended to support, contour, and provide structure for the chin. Made of high density polyethylene material with an interconnecting pore structure that supports tissue ingrowth. The device is available in various shapes and can be cut by the surgeon to meet the patient's specific needs.	78 x 9mm	Transferred to S2001
MW011	PMT Integra tissue expander	Tissue/skin Expander	5cc to 1400cc	Transferred to S2008

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
<b>Medtronic Australasia Pty Ltd</b>				
<b>Cardiac</b>				
MC812	Consulta CRT-D Model D234TRK	Fully Automatic, Wireless Implantable cardioverter defibrillator with cardiac resynchronization therapy (CRT), and therapies for ventricular and atrial tachyarrhythmia	38cc, 68g	Deleted
MC813	Secura VR Model D234VRG	Fully Automatic, Wireless Implantable cardioverter defibrillator with high energy capability and therapies for ventricular tachyarrhythmia	37cc, 68g	Deleted
MC814	Secura DR Model D234DRG	Fully Automatic, Wireless Implantable cardioverter defibrillator with high energy capability and therapies for ventricular and atrial tachyarrhythmia	37cc, 68g	Deleted
MC983	Consulta CRT-D Model D214TRM	Fully Automatic, Wireless Implantable cardioverter defibrillator with cardiac resynchronization therapy (CRT), and therapies for ventricular and atrial tachyarrhythmia. DF-4 lead connector.	38cc	Deleted
MI015	PROTECTA XT VR D354VRG	Single chamber implantable cardioverter defibrillator with SmartShock Technology, OptiVol 2.0 Fluid Status Monitoring and Complete Capture Management Specifications (DF-1)	Volume 37cc, Mass 68g	Deleted
MI016	PROTECTA XT DR D354DRG	Dual chamber implantable cardioverter defibrillator with SmartShock Technology, OptiVol 2.0 Fluid Status Monitoring and Complete Capture Management Specifications (DF-1)	Volume 37cc, Mass 68g	Deleted
MI017	Protecta XT CRT-D D354TRG	Implantable Cardioverter Defibrillator with Cardiac Resynchronisation therapy, SmartShock Technology, OptiVol 2.0 Fluid Status Monitoring and Complete Capture Management Specifications (DF-1)	Volume 38cc, Mass 68g	Deleted
MI024	Protecta XT VR D354VRM	Single chamber implantable cardioverter defibrillator with SmartShock Technology, OptiVol 2.0 Fluid Status Monitoring and Complete Capture Management Specifications (DF-4)	Volume 38cc, Mass 73g	Deleted
MI025	Protecta XT DR D354DRM	Dual chamber implantable cardioverter defibrillator with SmartShock Technology, OptiVol 2.0 Fluid Status Monitoring and Complete Capture Management Specifications (DF-4)	Volume 41cc, Mass 73g	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
MI026	Protecta XT CRT-D D354TRM	Implantable Cardioverter Defibrillator with Cardiac Resynchronisation therapy, SmartShock Technology, OptiVol 2.0 Fluid Status Monitoring and Complete Capture Management Specifications (DF-4)	Volume 40cc, Mass 73g	Deleted
MI034	Secura DR Model D214DRM	Fully Automatic, Wireless Implantable cardioverter defibrillator with high energy capability and therapies for ventricular and atrial tachyarrhythmia. DF-4 lead connector.	37cc	Deleted
MI054	Secura VR Model D214VRM	Fully Automatic, Wireless Implantable cardioverter defibrillator with high energy capability and therapies for ventricular tachyarrhythmia. DF-4 lead connector.	37cc	Deleted

## Morton Surgical Pty Ltd ta Matrix Surgical Company

### General Miscellaneous

MG002	Heliogast Adjustable Gastric Band	Silicone and Titanium	HAGE (evolution), small band - inner diameter 26mm; HAGA (advanced), large band - inner diameter 28mm	Deleted
-------	-----------------------------------	-----------------------	---	---------

## N Stenning & Co Pty Ltd

### General Miscellaneous

NG068	Embozene Tandem Microspheres	Embozene TANDEM Microspheres are spherical tightly calibrated, biocompatible, non-resorbable, hydrogel microspheres with an inorganic polyzene-F coating, colour coded to indicate microsphere size. Each 20mL syringe is filled with 2 ml of product suspended in a non-pyrogenic physiological saline transport solution. Embozene TANDEM Microspheres may be loaded with Doxorubican-HCl and Irinotecan-HCl up to 50mg/ml which can elute local, controlled, sustained doses of drugs to targeted tumour sites after embolisation.	Tandem 2ml 40 \u00B1 10 \u00B5m Tandem 2ml 75 \u00B1 15 \u00B5m Tandem 2ml 100 \u00B1 25 \u00B5m	Transferred to BS315
NG069	Embozene Tandem Microspheres	Embozene TANDEM Microspheres are spherical tightly calibrated, biocompatible, non-resorbable, hydrogel microspheres with an inorganic polyzene-F coating, colour coded to indicate microsphere size. Each 20 mL syringe is prefilled with 3 mL of product suspended in a non-pyrogenic physiological saline transport solution. Embozene TANDEM microspheres may be loaded with Doxorubican-HCl and Irinotecan-HCl up to 50 mg/ml which can elute local, controlled, sustained doses of the drugs to targeted tumour sites after embolisation.	Tandem 3ml 40 \u00B1 10 \u00B5m Tandem 3ml 75 \u00B1 15 \u00B5m Tandem 3ml 100 \u00B1 25 \u00B5m	Transferred to BS316

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
<b>Vascular</b>				
NG067	Embozene Microspheres	Uniquely colour-coded microspheres used in the embolization of hypervascular tumours and arteriovenous malformations	In between the range of 40 - 1130 micron - 2ml	Transferred to BS314
<b>Orthofix Australia Pty Limited</b>				
<b>Spinal</b>				
OF011	Blackstone Medical Inc Anterior Cervical Plate System	Bone Screw, Primary/Rescue Screw	12, 14, 16mm	Deleted
OF019	Blackstone Top Locking Plate	Holds bone screws and anterior cervical plate in place	One size only	Deleted
OF028	Blackstone Medical Inc Anterior Cervical Plate System	Titanium alloy plate L1-L2	22mm - 54mm	Deleted
OF030	Blackstone Medical Inc Anterior Cervical Plate System	Titanium alloy plate L3	56mm - 90mm	Deleted
OF040	Blackstone - Fixed System Washer	Titanium	One Size Only	Deleted
OF043	Blackstone - Staple	Titanium	2/4/6mm	Deleted
OF049	Blackstone Medical Pillar TL	Transforaminal Lumbar Interbody Fusion PEEK cage	9mm x 27mm, 0 degrees (parallel). 8mm -16mm height in 1mm increments. 9mm x 27mm, 8 degrees (lordotic). 8mm -16mm height in 1mm increments	Deleted
OF050	Blackstone Medical Pillar PL	Posterior Lumbar Interbody Fusion PEEK Cage	9mm x 23mm, 0 degrees (parallel). 8mm - 16mm height in 1mm increments. 9mm x 23mm, 8 degrees (lordotic). 8mm-16mm height in 1mm increments	Deleted
OF053	Blackstone - Spinal Fixation System - Mono Axial Screws	Composition - Titanium	5.5mm x 40mm, 5.5mm x 45mm, 6.5mm x 40mm and 6.5mm x 45mm	Deleted
<b>Orthotech Pty Ltd</b>				
<b>Hip</b>				
OH279	JRI Furlong Total Hip System	Uncemented Femoral Stem, titanium, grit blast, with HA,	12 /14 Taper, and 5mm offset, short and standard spigot	Transferred to TG024
OH280	JRI Furlong active stem	JRI Furlong Active Stem, Titanium Alloy with Supravit Coating, 12/14 Taper	12 / 14 Taper, 9-16mm	Transferred to TG008
OH281	JRI Furlong revision Total Hip System	Modular femoral stem, Grit blast, HA	12 / 14 Taper, 200mm long, 200ml long 5mm offset	Transferred to TG019
OH282	JRI Total Hip System	JRI Furlong Revision Stems, Titanium Alloy with Supravit Coating, 12 / 14 Taper	12 / 14 taper, 250mm long, 250ml long 5mm offset	Transferred to TG017
OH283	JRI Total Hip System	JRI Furlong Revision Plus Stem, Titanium Alloy with Supravit Zoned Coating, 12 / 14 Taper	12 / 14 Taper, 200mm long and 200ml long with 5mm offset	Transferred to TG013
OH284	JRI Total Hip System	JRI Furlong Revision Plus Stem, Titanium Alloy with Supravit Zoned Coating, 12 / 14 Taper	12 / 14 Taper, 250mm long, and 250ml long with 5ml offset	Transferred to TG002
OH286	JRI CSF Cup	Grit blast with HA, titanium	44mm - 70mm	Transferred to TG004

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
OH287	JRI CSF Cup insert	UHMWPE unmodified	28mm, 32mm, 36mm	Transferred to TG009
OH288	JRI CSF Cup Ceramic Insert	Ceramic Inserts, Alumina oxide, BioloX Forte	28mm	Transferred to TG012
OH289	JRI CSF Hip System	CSF Bone Screws, Titanium Alloy	15mm - 50mm	Transferred to TG001
OH290	JRI Cemented Cup	UHMWPE Cup unmodified Cemented	28mm and 32mm	Transferred to TG007
OH292	JRI Ceramic femoral Head	Total Hip Replacement (THR) Ceramic Femoral Heads, BioloX Delta Ceramic, and Alumina Oxide	12 / 14 Taper, BioloX Delta Ceramic head sizes 28mm, 32mm. Alumina Oxide head sizes 28mm	Transferred to TG020
OH293	JRI Femoral Head	Cobalt Chrome Femoral; Heads	28mm, and 32mm head sizes	Transferred to TG005
OH294	JRI Total Hip System	Ceramic femoral head with sleeve	28, 32, 36mm	Transferred to TG015
OH295	JRI CSF Plus Cup	Titanium, plasma, HA for ceramic insert	44mm - 70mm O/D	Transferred to TG027
OH296	JRI CSF Plus Cup	Titanium, plasma, HA for modified poly insert	44mm - 70mm O/D	Transferred to TG022
OH297	JRI CSF Plus Cup ceramic Insert	CSF Plus Cup Liner, Delta Liner	28mm, 32mm, 36mm	Transferred to TG011
OH298	JRI CSF Plus Cup Insert	UHMWPE modified	10 degree liners 28mm, 32mm, 36mm	Transferred to TG025
OH300	JRI Total Hip System	Total Hip Replacement (THR) Ceramic Femoral Heads, BioloX Delta Ceramic, and Alumina Oxide	12 / 14 Taper, BioloX Delta Ceramic head sizes 36mm & 40mm, Alumina Oxide head sizes 28mm	Transferred to TG026
OH378	Securus Modular Revision Hip System	Curved Femoral Stem	250 - 300mm	Transferred to TG023
OH379	Securus Modular Revision Hip System	Modular Proximal Body	Std & 5mm offset	Transferred to TG018
OH380	Securus Modular Revision Hip System	Modular Body Cap	One size	Transferred to TG016
OH381	Securus Modular Revision Hip System	Modular Locking Bolt	One Size	Transferred to TG006
OH382	Securus Modular Revision Hip System	Distal locking screws	30-70mm	Transferred to TG010
OH397	Furlong Hemi Range	Cemented femoral stem	XXS, XS, S, M & L	Transferred to TG014
OH398	Furlong Hemi Range	Bipolar head	40-58mm	Transferred to TG003
OH520	Furlong Evolution	Collarless and collared 126-133 (degree)	6-17	Transferred to TG021

### Specialist Orthopaedic

OH488	DynafitSystem Staples	Foot Staple	8-10mm 90o and 26o	Expanded into OH527, OH528
OH504	G21 Bone Cement with Antibiotic+ HiVac7 delivery system	Bone cement with antibiotic \u2013 low & std viscosity	20-60g	Expanded into OH529, OH530
OH505	G21 Bone Cement + HiVac7 delivery system	Bone cement \u2013 low & std viscosity	20-60g	Expanded into OH525, OH526

### Pyramed Pty Limited

#### Cardiac

PY034	Figulla Flex II PFO Occluder	The Figulla Flex II PFO is an implant for closure of a patent foramen ovale (PFO) defect in the heart and is deployed via a minimally invasive catheter delivery system. The Figulla Flex PFO is intended for patients with substantiated Right-Left-Shunt in the contrast agent echocardiography under Valsalva maneuver and presence of a "cryptogenic", paradoxical thromboembolism.	Figulla Flex II PFO Disc (mm): 23/25, 16/18, 27/30, 31/35	Transferred to ME197
-------	------------------------------	---	---	----------------------

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
PY035	Figulla Flex II ASD Occluder	Occlutech Flexx II ASD Occlusion device	Figulla Flex II ASD Waist \u00D8 (mm): Size 4,5,6,7,5,9,10,5,12,13,5,15,16,5,18,21,24,27,30,3 3,36,39,40	Transferred to ME196

## RQ Solutions Medical Devices Distribution Support

### Specialist Orthopaedic

RQ001	Orthomed Metallic Snap-off screw	The Snap-off screw is indicated for fixation of bone fragments. Snap off screws are breakable self tapping screws, which present a breakable cartridge necessary to the screw setting. They are particularly intended to repair the 2nd, 3rd, 4th and 5th metatarsal of the foot.	Thread diameter: 2mm/cortical thread Core diameter: 1mm Overall length: depending on reference Threaded length: depending on reference 2 opposite slots at the tip (self drilling and self tapering profile) 4 slots in the screw head for screwdriver grip	Transferred to ER298
RQ002	Orthomed Interference Screws with Atraumatic Thread	The screws are intended to fix an autogenous transplant (hamstring tendons) or a synthetic graft in the femur and tibia during an Anterior Cruciate Ligament (ACL) reconstruction.	9 sizes Diameter: 7,8 and 9mm Length: 18, 23 and 28mm	Transferred to ER300
RQ003	Orthomed Metallic Compression Screw	The compression screw is indicated for fixation of bone fragments for forefoot surgery. The compression screws are self-tapping cannulated screws, which present two types of thread diameter (one on the screw head and a body part of the screw). They are particularly intended to repair the 1st metatarsal in the foot. The pitch difference between the proximal and the distal part allows creating a compression are between the upper and lower bone fragment of 1st metatarsal improving the healing process.	10 sizes Head diameter: 3.8mm Proximal diameter: 3mm Internal diameter: 1.10mm Different length from 10mm to 34mm by pitch of 2mm between each reference	Transferred to ER294
RQ006	Orthomed MISBIO Bioresorbable Interference Screws	The Orthomed Misbio Bioresorbable interference screws are intended to set an autogenous transplant (hamstring tendons or patellar tendon) in the femur and tibia during an Anterior Cruciate Ligament (ACL) reconstruction	Diameter 6-10mm, length 23mm Diameter 7-10mm, length 30mm Diameter 9-11mm, length 35mm Details: MISBIO 6.23: 6mm diameter, 23mm length MISBIO 7.23: 7mm diameter, 23mm length MISBIO 8.23: 8mm diameter, 23mm length MISBIO 9.23: 9mm diameter, 23mm length MISBIO 10.23: 10mm diameter, 23mm length MISBIO 7.30: 7mm diameter, 30mm length MISBIO 8.30: 8mm diameter, 30mm length MISBIO 9.30: 9mm diameter, 30mm length MISBIO 10.30: 10mm diameter, 30mm length MISBIO 9.35: 9 mm diameter, 35mm length MISBIO 10.35: 10mm diameter, 35mm length MISBIO 11.35: 11mm diameter, 35mm length	Transferred to ER303
RQ007	Orthomed Interference Screw Bone to Bone	The screws are intended to set an autogenous transplant (hamstring tendons or patellar tendons) in femur and tibia during an Anterior Cruciate	9 sizes: Diameter: 7, 8 & 9mm Lengths: 18, 23 & 28mm	Transferred to ER299

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
RQ070	LIGASTIC reinforcement ligament	<p>The LIGASTIC reinforcement artificial ligament is to constitute a reinforcement or support during healing &amp; repair of the Patellar tendon, achilles tendon, antero lateral laxity of ankle, medial collateral ligament and posterior cruciate ligament. There are different products for the reinforcement of the collateral ligaments, anterior cruciate ligament, achilles tendon, shoulder, patella tendon, popliteal tendon, ankle anterior-lateral laxity and acromial-clavicular dislocation.</p> <p>The product consists of warp-knitted fibres of PET. It is knitted in bands of 8, 15, 18, 22, 36 or 72 warp threads.</p> <p>The basic structure of LIGASTIC is always identical. However, the presentation varies depending on the actual joint affected. The ligaments presented are:</p> <ul style="list-style-type: none"> <li>- cylindrical string, with variable length &amp; diameter. Diameter varies according to fiber content &amp; required resistance. Both ends of the ligament are soft, formed as a reinforced loop, allowing an easy passage through the bone tunnels. These ends have a temporary function &amp; are able to be cut off after the tunnel passage.</li> <li>- band, for reinforcement of a flat anatomical structure such as the medial collateral ligament in the knee</li> <li>- reinforcement ligament (uncoated LIGASTIC) - for support and healing of intra-articular or extra-articular natural ligaments</li> </ul>	<p>7 sizes:</p> <p>LCA60NEF: 30 threads, length 300mm, diameter 6mm</p> <p>LT60:30 threads, length 290mm, width 6 or 8.8mm</p> <p>RTF22: 11 threads, length 360mm, width 8mm</p> <p>RTF18: 9 threads, length 360mm, width 7mm</p> <p>LCP88NE: 22 threads, length 550mm, width 5.5mm</p> <p>LLI30: 15 threads, length 350mm, width 22.5mm</p> <p>L30: 15 threads, length 360mm, diameter 4.5mm</p> <p>LT60: 30 threads, length 290mm, width 8.8mm</p>	Transferred to ER301

## Signature Orthopaedics

### Hip

SF040	Global Acetabular Cancellous Screw	Titanium 6.5mm cancellous bone screw	15mm to 50mm	Deleted
-------	------------------------------------	--------------------------------------	--------------	---------

### Spinal

SF002	VerteLoc\u2122 Spine System	Interbody Fusion Device	7mm - 18mm	Deleted
SF004	ZygoLoc	Facet Fusion Cage	3-15mm	Deleted
SF010	Centaur Screw System	Set Screw	small, standard and pelvic	Deleted
SF013	Centaur Lateral Connector, Centaur Screw System	Cross Connector - Rod to Rod	20-500mm	Deleted
SF020	Aurora SPINE System	Thoracolumbar plate	15mm-51mm	Deleted
SF051	Aquila Plate System	Anterior Cervical Plate	1 level and 2 level plates (21 - 55mm Length)	Deleted
SF052	Aquila Plate System	Anterior Cervical Plate	2, 3 and 4 level cervical plates (55m - 109mm length)	Deleted
SF065	Aquila Plate System	Bone Screw	10mm-20mm	Deleted
SF066	Aquila Plate System	Bone Screw	10mm-20mm	Deleted



Billing Code	Product Name	Further Description	Size	Explanation of Deletion
SF071	Centaur Screw System	Standard Rod	30mm-480mm	Deleted
SF073	VerteLoc Spine System	Cervical Cage IF	Extra small - Extra Large	Deleted
SF074	VerteLoc Spine System	Bone Screw	10-20mm	Deleted
SF075	VerteLoc Spine System	Bone Screw	10-20mm	Deleted
SF080	Centaur Screw System	Polyaxial Screw	4mm - 13mm	Deleted
SF085	Australis Spine System	Thoracolumbar Cage	Small to extra large	Deleted

## Smith & Nephew Pty Ltd

### Knee

SN836	Journey Knee System Tibial Uni Baseplate	Uni Tibial Baseplate,Titanium, cemented, for fixed insert	1-8	Deleted
-------	--	---	-----	---------

### Specialist Orthopaedic

SL024	RCI PK	Interference screw	5x22, 5x22L, 6x22, 6x22L, 6x27, 6x27L, 7x22, 7x22L, 7x27, 7x27L, 7x32, 7x32L, 8x22, 8x22L, 8x27, 8x27L, 8x32, 9x27, 9x27L, 9x32, 9x37, 10x32, 10x37, 11x32, 11x37, 12x32, 12x37	Deleted
-------	--------	--------------------	---	---------

## St Jude Medical Australia Pty Ltd

### Cardiothoracic

SJ316	Seguin Annuloplasty Ring	Semi-rigid annuloplasty ring	24-40mm	Deleted
-------	--------------------------	------------------------------	---------	---------

### Neurosurgical

SJ133	Eon Mini Neurostimulation System - Implantable Pulse Generator 3788	Eon Mini Neurostimulation System - Implantable Pulse Generator 3788	18cc, 29g	Deleted
SJ182	Eon Neurostimulation System	Eon Neurostimulation Implantable Pulse Generator #3716	Height = 59mm	Deleted
SJ203	Brio DBS Charging System 6721	Charging System Model 6721	One size only	Deleted
SJ211	Eon Neurostimulation System - Mobile Charging System 3711	Mobile Charging System	One size	Deleted
SJ347	Prodigy Neurostimulator	Prodigy Neurostimulator	29g	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
--------------	--------------	---------------------	------	-------------------------

## Teleflex Medical Australia

### Ear, Nose & Throat

TX048	Ventilation Tubes / Grommet	a) Tubingen Gold Type b) Tubingen Titanium Type c) T Tube d) Reuter Bobbin Type e) Paparella Type f) Collar Button Type g) Star Tube h) Bevel Button Type i) Straight Tube j) Armstrong Type k) Shepard Type l) Donaldson Type m) Angular Type n) Myringotomy Tube with Trocar o) Ventilation Tubes with Eyelets p) Minimal Type	a) I.D. 1.25mm w/w, I.D.1.25mm wo/w, I.D. 1.5mm wo/w b) one size only c) L- 6, 7.5, 12 d) I.F.O.D 2.5mm, 2.7mm e) one size only f) one side only g) one size only h) one size only i) i.d.-0.9mm,1.1mm,1.15mm j) L-3.8mm, 3.7mm, 10mm k) L- 2.2mm w/t, 2.2mm wo/t,2.4mm w/t l) S- one size only, F one size only m) one size only n) Silver or Titanium o) one or two p) one size only	Transferred to A2004
TX049	Tracheal Speaking Valve	1. Passy -Muir - a) Purple b) Clear c) Aqua 2. Provox - a) Two b) Vega c) Vega Puncture Set	1. One size only 2. a) 4.5, 6, 8, 10, 12.5, 15mm b) 17, 20, 22.5Fr in 4, 6, 8, 10, 12.5, 15mm c) 17, 20, 22.5Fr in 8, 10, 12.5, 15mm	Transferred to A2002
TX050	Vocal Chord Medialisation Implant	Montgomery Laryngeal Keel	12, 14, 16 - Clear and Radiopaque	Transferred to A2001
TX051	Thyroplasty	Montgomery Thyroplasty Implant System	male - 8,9,10,11,12,13 female - 6,7,8,9,10	Transferred to A2005
TX052	Extra-tracheal Cannula	a) Montgomery Short-Term Cannula b) Montgomery Long-Term Cannula	a) 4, 6, 8, 10 b) 4, 6, 8, 10 for Stoma Sizes 17-25, 25-33, 33-41, 41-49, 49-57, 57-65mm	Transferred to A2006
TX053	Intra-tracheal Cannula	Montgomery Safe-T-Tubes Clear and Radiopaque except c) a) Paediatric b) Standard c) Hebeler d) Thoracic e) Extra-Long f) Tapered g) Montgomery Salivary Bypass Tube h) Har-el Pharyngeal Tube i) Montgomery Esophageal Tube j) Trachael Fenestrator k) Moore Trache	a) 6, 7, 8, 9 b) 10, 11, 12, 13, 14, 15, 16 c) 12, 15, 18 d) 10, 11, 12, 13, 14, 15, 16 e) 10, 11, 12, 13, 14, 15, 16 f) 8/10, 10/13 g) 8, 10, 12, 14, 16, 18, 20 h) one size only i) one size only j) size 6, 8 and 10 k) size 6 and 8	Transferred to A2003

### General Miscellaneous

TX038	EndoMaxx Fully Covered Oesophageal Stent	Oesophageal stent, radiopaque, self expanding nitinol covered and premounted on delivery system	Diameter 19mm-23mm. Length 70mm-150mm.	Deleted
TX039	Alimaxx-B Uncovered Biliary Stent System	Biliary stent, radiopaque, self expanding nitinol and premounted on delivery system	Diameter 8mm-10mm. Length 40mm-80mm	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
TX041	ALIMAXX-E Esophageal Stent System	Esophageal stent, radiopaque, self expanding nitinol covered stent, premounted on a delivery system	Diameter 18-22mm, Length 70-120mm	Deleted
TX042	AERO Tracheobronchial Stent System	Radiopaque self-expanding nitinol stent and delivery system	Diameter 10-20mm, length 20-80mm	Deleted

### Urogenital

TX044	ReSolve (Nephrostomy) Locking Drainage Catheter	Nephrostomy Catheter	6.5F x 15cm, 7.5F x 15cm, 8.5F x 15cm, 10F x 15cm, 12F X 15cm, 14F x 15cm, 6.5F x 25cm, 7.5F x 25cm, 8.5F x 25cm, 10F x 25cm, 12F x 25cm, 14F x 25cm, 6.5F x 40cm, 7.5F x 40cm, 8.5F x 40cm, 10F x 40cm, 12F x 40cm, 14F x 40cm	Deleted
-------	---	----------------------	---	---------

### Terumo Australia

#### Vascular

TU039	Rapidax\u2122 Rapid Access Haemodialysis Graft	Triple layered self sealing haemodialysis Access Graft. ePTFE internal and external layer with proprietary elastomer inner layer - Straight and Pre-curved	6mm x 20-40cm	Deleted
TU051	AZUR Peripheral HydroCoil, pushable	AZUR Peripheral HydroCoil, pushable	2-20mm diameter, length 2-15cm	Deleted

### W L Gore & Associates

#### General Miscellaneous

GT266	Gore\u00AEBio-A\u00AE Fistula Plug	Gore\u00AEBio-A\u00AE Fistula Plug	Six 9cm Tubes, 16mm disk	Deleted
-------	------------------------------------	------------------------------------	--------------------------	---------

### WILLIAM A. COOK AUSTRALIA PTY. LTD.

#### General Miscellaneous

WC174	Zilver Biliary Stent System	Self expanding Metal Biliary Stent System	Diameters: 6, 8, 10mm; Lengths: 4, 6,8cm	Deleted
-------	-----------------------------	---	--	---------

### Urogenital

WC027	Polyurethane Double Pigtail Stent Set	Polyurethane Double Pigtail Stent Set	Multiple French Sizes and Lengths	Deleted
WC131	Endo-sof double pigtail stent Set with AQ Hydrophilic Coating and Sideholes	Endo-sof double pigtail stent Set with AQ Hydrophilic Coating and Sideholes	Multiple French Sizes and Lengths	Deleted
WC137	Pediatric Urethral C-Stent	Pediatric Urethral C-Stent	Multiple French Sizes, Length 15cm	Deleted

### Vascular

WC301	Zenith Low Profile AAA Endovascular Graft - Contralateral/Ipsilateral Iliac Leg	Contralateral/Ipsilateral Iliac Leg	Diameter: 10-24mm, Working Length: 36-120mm	Deleted
-------	---	-------------------------------------	---	---------

### Wright Medical Australia

#### Hip

WR084	Femoral Head	Femoral Head Cast Cobalt Chrome Alloy	22.5mm - 32mm Femoral Head and neck ranging from -3.5mm to +7mm	Transferred to M1035
-------	--------------	---------------------------------------	---	----------------------

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
WR085	Profemur Xm Femoral Stem	Femoral Stem, CoCr, Cemented, Polished	0, 1, 2, 3, 4	Transferred to M1034
WR095	Procotyl L Acetabular Cup	Beaded, Titanium HA coated cup	42 - 68mm	Transferred to M1033
WR099	Procotyl L Delta Ceramic Liner	Delta ceramic liners	28, 32, 36, 40, 44 ID	Transferred to M1029
WR100	Delta Modular Heads	Biolox Delta Heads	36mm - 44mm in -3.5mm, 0mm, +3.5 neck lengths	Transferred to M1028
WR101	Delta Modular Heads	Biolox Delta Heads	28mm, 32mm in -3.5mm, 0mm, + 3.5 neck lengths	Transferred to M1038
WR103	Dynasty Acetabular System /poly liner	A-class Crosslinked polyethylene liner, standard and lipped, large head options	ID 22,28,32,36, 38,40,42,46	Transferred to M1025
WR106	Procotyl L A-class Crosslinked polyethylene liner,standard & lipped, large head options	A-class crosslinked poly liners, standard and lipped liners	28, 32, 36, 40 ID	Transferred to M1026
WR114	Dynasty Acetabular System/ Actebular Shell	Beaded Acetabular shell non HA	42 - 76mm in 2mm increments	Transferred to M1032
WR118	WMT Perfecta Hip System	PMMA cement Plug	9 - 21mm in 1mm increments	Transferred to M1030
WR139	Profemur Modular Neck	Titanium and CoCr Femoral Necks	Length 27 - 38.5mm; Neutral to Ante/Retro - Varus Valgus 2	Transferred to M1041
WR222	Profemur L	Profemur L Hip Stem	Size 1 - 11	Transferred to M1031
WR249	Cement Cetraliser	Profemur XM Cement Centraliser	N/A	Transferred to M1039
WR253	Profemur Hip System	Primary stem, tapered wedge, plasma coated	1-12	Transferred to M1024
WR254	Profemur Hip System	Primary stem, tapered wedge, plasma coated, fixed neck	1-12	Transferred to M1040
WR265	Transcen/Lineagte Articulation System	Metal Femoral Head CoCr	28mm- 36mm, 4 lengths	Transferred to M1027
WR271	Bone Screw	Cancellous Bone Screw Titanium	6.5mm x 1.5 \u2013 6 cm in 0.5cm increments	Transferred to M1036
<b>Knee</b>				
WR161	Medial Pivot Knee System	Stemmed Condylar medial pivot component, uncemented, beaded, non HA, CoCr used with anterior lipped insert, minimally constrained	1, 2, 3, 4, 5, All in left & right	Deleted
WR168	Medial Pivot Knee System	Revision tibial component, titanium, cemented with offset stem	1, 1+, 2, 2+, 3, 3+, 4, 4+, 5, 5+ in Left & Right	Deleted
WR169	Medial Pivot Knee System	Tibial component, beaded, titanium	1, 1+, 2, 2+, 3, 3+, 4, 4+, 5, 5+, 6	Deleted
WR178	Medial Pivot Knee System	Offset tibial block augments, screw fixation for use with offset stemmed tibial component	1, 2/1, 3/2, 4/3, 5/4, 6/5, All in 5mm, 10mm, 15mm	Deleted
WR180	Medial Pivot Knee System	Stemmed condular medial pivot component, cemented, grit blasted, minimally stabilised, femoral component, CoCr, usually used with a anterior lipped insert	1, 2, 3, 4, 5 All in Left & Right	Deleted
WR181	Medial Pivot Knee System	Tibial component, beaded, titanium, HA	1, 1+, 2, 2+, 3, 3+, 4, 4+, 5, 5+, 6	Deleted
WR183	Medial Pivot Knee System	Uncemented femoral component, beaded, HA, CoCr, usually used with a anterior lipped insert	1 - 6 Left and Right	Deleted
WR184	Medial pivot KneeSystem	Uncemented, beaded, non HA, CoCr, usually used with a anterior lipped insert	1, 2, 3, 4, 5, 6	Deleted

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
WR185	Medial pivot Knee System	Cemented, grit blasted, totally constrained, femoral component, CoCr	1, 2, 3, 4, 5	Deleted
WR187	Medial Pivot Knee System	totally constrained, UHMWPE, fixed tibial insert with reinforcing peg	1, 2, 3, 4, 5, 6 in 10 - 30mm in 2mm increments	Deleted
WR188	Medial Pivot Knee System	Tibial component, CoCr, cemented	1, 1+, 2, 2+, 3, 3+, 4, 4+, 5, 5+, 6, 7, 8, 8+	Deleted
WR189	Medial Pivot Knee System	Posterior stabilised, fixed tibial insert, UHMWPE	1, 2, 3, 4, 5, 6, 7, 8, 10-24mm	Deleted
WR190	Medial pivot Knee System	All poly patella, recessed and onlay, single peg or tri-peg	25, 28, 32, 35, 38, 41mm	Deleted
WR193	Medial Pivot Knee System	Femoral Augments titanium, separate distal and posterior with screw attachment	1 - 5, 5mm thick, 10mm thick, posterior and distal	Deleted
WR195	Medial Pivot Knee System	Uncemented stems, femoral and tibial	11-22mm, 100-140mm lengths	Deleted
WR197	Medial Pivot Knee System	Tibial augments, wedges and blocks, screw attachment	15 degree, sizes 1, 2/1+, 3/2+, 4/3+, 5/4+, 6/5+, Blocks 5mm, 10mm, 15mm	Deleted
WR198	Medial Pivot Knee System	Revision tibial component, titanium, cemented	0+, 1, 1+, 2, 2+, 3, 3+, 4, 4+, 5, 5+, 6	Deleted
WR199	Medial Pivot Knee System	Cemented stems, femoral and tibial	10, 12, 14, 16, 18 x 65mm, 15 x 30mm	Deleted
WR200	Medial Pivot Knee System	Option bone screw for tibial component	6.5mm x 1.5 - 6cm in 0.5cm increments long	Deleted
WR223	Medial Pivot Knee System	Minimally stabilised, fixed, tibial insert, UHMWPE	1, 1+, 2, 2+, 3, 3+, 4, 4+, 5, 5+, 6, 6+, 7, 7+, 8, in 10 - 24 mm thicknesses	Deleted
WR251	Advance Medial Pivot Stem Cap	Modular Tibial Base Stem Cap	N/A	Transferred to M1037
WR262	Medial pivot Knee System	Cemented, grit blasted, minimally stabilised, femoral component	1, 2, 3, 4, 5, 6, 7 & 8	Deleted
WR263	Medial pivot Knee System	Cemented, grit blasted, prosterior stabilised, femoral component	1, 2, 3, 4, 5, 6, 7 & 8	Deleted

## Zimmer Biomet Pty Ltd

### Hip

BI053	Bi-Metric Hip System	Revision 250 Stem, Titanium, plasma, non HA	9-17mm	Deleted
BI054	Bi-Metric Hip System	Revision 300mm Femoral Stem, Titanium, plasma non HA	9-17mm	Deleted
BI110	Mallory Head Hip System	Acetabular Shell	40-80mm	Compressed into ZI939
BI112	Mallory Head Hip System	Acetabular shell, solid, plasma, HA	40-80mm	Compressed into ZI939
BI896	M2A Magnum System	Femoral Head, Cobalt Chrome, metal on metal	38-60mm	Deleted
ZI055	Ceramic Femoral Head	Ceramic Head 12/14 Taper, Biolox Forte	28 - 36mm	Deleted
ZI311	Spacer-G	Temporary implantable device for the treatment of infected total hip replacement, Composition - stainless steel core surrounded by PMMA impregnated with Gentamicin	Short Stems, Head size 46mm, 54mm, 60mm	Transferred to ER295
ZI312	Spacer-G	Temporary implantable device for the treatment of infected total hip replacement, Composition - stainless steel core surrounded by PMMA impregnated with Gentamicin	Long Stems, Head size 46mm, 54mm, 60mm	Transferred to ER296

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
ZI381	Ceramic Head - Sulox	Head - Ceramic	22 to 32mm	Deleted
<b>Knee</b>				
BI945	Vanguard PFR	Femoral Component, Ti	XS to Xlarge	Deleted
ZI313	Spacer-K	Temporary implantable device for the treatment of infected total knee replacement, Composition - stainless steel core surrounded by PMMA impregnated with Gentamicin	Small, Medium, Large	Transferred to ER297
<b>Plastic and Reconstructive</b>				
BI448	Lorenz 2.4mm Locking Reconstruction System	Plate, curved, std or low profile, 2.4mm system, titanium	16-32 hole	Compressed into ZI938
BI452	Lorenz 2.4mm Locking Reconstruction System	Plate, curved. Std or low profile, 2.4mm System, Titanium	16 hole	Compressed into ZI938
BI455	Lorenz 2.4mm Locking Reconstruction System	Plate, Double Angled, 2.0 or 2.4mm System, Titanium	5 x (22 or 24) x 5 hole	Compressed into ZI938
BI986	Lorenz 2.4mm Locking Reconstruction System	Plates Plate, Angled 2.4mm System Ti	5 x (22 -24) hole	Compressed into ZI938
<b>Specialist Orthopaedic</b>				
BH112	Optipac	Refobacin Plus, Prefilled Vaccum Mixing System, Gentamicin Sulphate, PMMA Cement	40gm	Compressed into ZI937
BH113	Optipac	Refobacin Plus, Prefilled Vacuum mixing system, Gentamicin Sulphate, PMMA Cement	60-80gm	Compressed into ZI937
BH239	Versa Nail	TTC Fusion Nail	10mm x 15cm up to 12mm x 30cm	Deleted
BH248	Compression Hip Screw System	Stainless Steel Trochanteric Support Plate	Small-Large	Deleted
BH249	Compression Hip Screw System	Stainless Steel Trochanteric Support Plate	Small-Large	Deleted
BH250	Compression Hip Screw System	Stainless Steel Trochanteric Support Plate	Small-Large	Deleted
BH267	VersaNail IM Nailing Platform	Tibial Nails, Titanium, Locking, cannulated.	8-13mm x 25 - 50cm	Deleted
BH304	Compression Hip Screw System	Compression hip screw plate, stainless steel	2 to 6 holes	Deleted
BH305	Compression Hip Screw System	Compression hip screw plate, stainless steel	7 to 15 holes	Deleted
BH306	Compression Hip Screw System	Compression hip screw plate, stainless steel	16 holes	Deleted
BH328	Versa Nail	TTC Nail Cap Unsterile Impinging and Non Impinging	One Size	Deleted
BH329	Versa Nail	TTC Nail Cap Sterile Impinging and Non Impinging	One Size	Deleted
ZI873	Cemex\u00AE System Genta	Bone Cement with Gentamicin packaged in a complex delivery system	40, 60, 70 and 80 grams	Transferred to ER304
ZI877	Cemex-genta bone cement	1 mix Methylmethacrylate/ Polymethylmethacrylate with gentamicin	1 mix 40g powder and 13.3g liquid	Transferred to ER302

**Grand Total: 459 Records**