

Prostheses List

Product Category	Assessment Body	Product Group	Product Sub-Group	Suffix	Billing Code	Min Benefit	Max Benefit	Sponsor	Product Name	Description	Size	ARTG	Last Updated	Note	Change1
02 - Ear, Nose & Throat	02.01 - EAR	02.01.05 - Ossicle/Middle Ear Prosthesis	02.01.05.05 - Active Middle Ear Implants	Coupler	US019	350		Med-El Implant Systems Australasia Pty Ltd	Vibrant Soundbridge System - Vibroplasty Couplers	Middle ear prostheses used with the Vibrant Soundbridge System VORP implant. Facilitate coupling between the FMT and a vibratory structure of the middle ear.	Functional length: <input type="checkbox"/> Vibroplasty-OW-Coupler - 2.50, 3.00, 3.50, 4.00mm <input type="checkbox"/> Vibroplasty-Clip-Coupler - 0.20, 0.50, 1.00, 1.50mm <input type="checkbox"/> Vibroplasty-Bell-Coupler - 0.20, 0.50, 1.00, 1.50mm <input type="checkbox"/> Vibroplasty-RW-Coupler - 0.85mm <input type="checkbox"/> Incus-Symphoni x-Coupler - Left and Right (2.5	188553	February 2017		New Listing
02 - Ear, Nose & Throat	02.01 - EAR	02.01.05 - Ossicle/Middle Ear Prosthesis	02.01.05.05 - Active Middle Ear Implants	Program mable, wireless - enabled	US020	6825		Med-El Implant Systems Australasia Pty Ltd	Vibrant Soundbridge System - Samba Audio Processor	Middle ear implant audio processor. The Samba audio processor is the externally worn component of the Vibrant Soundbridge System and connects to the implant through magnetic attraction. Fitting the Samba audio processor activates the Soundbridge.	Samba Lo, Samba Hi left and right versions, length = 35mm, width = 30mm, height = 10mm, weight <9g (inc battery)	270912	February 2017		New Listing
02 - Ear, Nose & Throat	02.01 - EAR	02.01.05 - Ossicle/Middle Ear Prosthesis	02.01.05.05 - Active Middle Ear Implants		US021	7470		Med-El Implant Systems Australasia Pty Ltd	Vibrant Soundbridge System - Vibrating Ossicular Prosthesis VORP503	Part of the Vibrant Soundbridge System, VORP503 is an active implant, which is implanted into the temporal bone under the skin. The VORP 503 consists of the Floating Mass Transducer (FMT), a conductor link, the electronics (demodulator), fixation wings and a magnet surrounded by a receiver coil.	One size (left and right) Length: 118mm, Width: 28mm, Thickness: 4.4mm, Weight: 10.1g	287904	February 2017		New Listing
10 - Vascular	10.01 - Vascular Stents	10.01.01 - Bare Metal Stents	10.01.01.03 - Self Expandable, ≥150mm		DL002	2970		DIVERSE DEVICES PTY LTD	iVolution ≥150mm	Self-expanding nitinol stent system	Diameters: 5-8mm <input type="checkbox"/> Lengths: 150-200mm <input type="checkbox"/> Catheter Length: 80-140cm	281201	February 2017		New Listing

Prostheses List

03 - General Miscellaneous	03.02 - DRUG DELIVERY DEVICES	03.02.02 - Infusion Pumps, Balloon Based	03.02.02.01 - Fixed Flow Rate	Set	CU003	260		Cosuda Pty Ltd	Neofuser Pump	Continuous flow rate ambulatory and disposable silicone balloon infusion pump	Volume 100ml, 275ml, 450ml	286946	February 2017		New Listing
03 - General Miscellaneous	03.02 - DRUG DELIVERY DEVICES	03.02.02 - Infusion Pumps, Balloon Based	03.02.02.02 - Variable Flow Rate	Set	CU004	375		Cosuda Pty Ltd	Neofuser Vario Pump	Variable flow rate ambulatory and disposable silicone balloon infusion pump	Volume 60ml, 100ml, 180ml, 275ml, 450ml	286946	February 2017		New Listing
03 - General Miscellaneous	03.02 - DRUG DELIVERY DEVICES	03.02.02 - Infusion Pumps, Balloon Based	03.02.02.02 - Variable Flow Rate	Bolus	CU005	450		Cosuda Pty Ltd	Neofuser Bolus Pump	Ambulatory and disposable silicone balloon infusion pump with bolus function	Volume 100ml, 275ml, 450ml	286946	February 2017		New Listing
02 - Ear, Nose & Throat	02.01 - EAR	02.01.05 - Ossicle/Middle Ear Prosthesis	02.01.05.05 - Active Middle Ear Implants		US018	7470		Med-El Implant Systems Australasia Pty Ltd	Vibrant Soundbridge System - Vibrating Ossicular Prosthesis VORP502x	An active direct-drive implantable middle ear hearing prosthesis, the VORP 502x consists of the Floating Mass Transducer (FMT), a conductor link, an electronic package and a magnet surrounded by an internal coil.	VORP502 A (Right), VORP502 B (Left) Length: 129mm, Width: 29mm, Thickness: 4.5mm, Weight: 9.3g	161702	February 2017		New Listing
End of work sheet															