

# Deleted Items for Prostheses List

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
<b>Device Technologies Australia Pty Ltd</b>				
<b>Ophthalmic</b>				
DE388	Oculentis LENTIS Mplus Toric LU-313MFT	Patient specific one-piece multifocal-toric plate haptic acrylic IOL	Overall Length 11.0mm and Optic Size 6.0mm + 0.0D to +36.0D inclusive in 0.01D increments. Cyl +0.25 to +6.99D inclusive in 0.01D increments	Expanded into DE494, DE495
DE392	Oculentis LENTIS LU-313-TY	Patient Specific: Foldable one-piece monofocal toric acrylic IOL. Violet Light Filter	Optic Size: 6.0mm, overall length 11.00mm. Sph +/-0.0D to +35.0D (in 0.01D increments) Cyl +0.25D to +6.99D cyl (in 0.01D increments)	Expanded into DE488, DE489
DE396	Oculentis LENTIS Mplus LS-313 MF	Foldable one-piece multifocal plate haptic acrylic IOL.	Overall Length 11.0mm and Optic Size 6.0mm Available diopters +-0.0D to +36.0D (0.05D steps)	Expanded into DE496, DE497
DE399	Oculentis LENTIS L-312	Foldable one-piece monofocal acrylic C Loop IOL	Optic size 6.0mm, Overall Length 12.00mm. Central Thickness 1.02mm (+22.0D). Available Diopters: -10.0D to 35.0D (1.0D steps). +10.5D to 29.5D (0.5D steps)	Expanded into DE490, DE491
DE400	Oculentis LENTIS L-313	Foldable one-piece monofocal acrylic plate haptic IOL	Optic size 6.0mm, Overall Length 11.00mm. Available Diopters: +-0.0D to 35.0D (1.0D steps). +10.5D to 29.5D (0.5D steps)	Expanded into DE492, DE493
DE401	Oculentis LENTIS L-402	Foldable three piece monofocal acrylic IOL.	Optic Size 6.0mm, Overall Length 13.00mm. +-0.0D to +35.0D (1.0D) + 15.5D to +26.5D (0.5D). Minus range on demand.	Expanded into DE484, DE485
DE402	Oculentis LENTIS LS-312Y	Foldable one-piece monofocal C Loop Acrylic IOL. Violet Light filter.	Optic Size 6.0mm, Overall Length 12.0mm. Available Diopters: +-0.0D to 35.0D (1.0D steps) +15.5D to 26.5D (0.5D steps).	Expanded into DE486, DE487
DE403	Oculentis LENTIS LS-313Y	Foldable one-piece monofocal acrylic plate haptic IOL. Violet Light Filter	Optic Size 6.0mm, Overall Length 11.00mm. Available Diopters: +-0.0D to 35.0D (1.0D steps). +15.5D to 26.5D (0.5D steps).	Expanded into DE498, DE499

## The trustee for Devitt Family Trust

### General Miscellaneous

TD001	Purple Surgical Ultimate circular stapler	The Purple Surgical Ultimate circular staplers can be applied throughout the alimentary tract to create end-to-end, end-to-side, and side-to-side anastomosis. The Purple Surgical Ultimate circular stapler places a staggered circular, double row of titanium staples in the tissue and resects the excess tissue, to create a circular anastomosis. The instrument can be activated when the handle is squeezed firmly as far as it can go. There are 5 diameters available, 24mm, 26mm, 29mm, 32mm, and 34mm with colour coded anvil heads.	24mm, 26mm, 29mm, 32mm, and 34mm.	Transferred to ME202
-------	---	--	-----------------------------------	----------------------

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
TD002	Purple Surgical Ultimate endoscopic linear cutter stapler reloading units, articulating	The Purple Surgical Ultimate endoscopic linear cutter stapler reloading units, articulating, are for single patient use only, and are used in conjunction with the Purple Surgical Ultimate endoscopic linear cutter stapler universal handles. A single use reloading unit can be replaced up to 24 times on the same universal handle, during the same procedure. The endoscopic linear cutter stapler articulating reloading unit is compatible with the endoscopic linear cutter stapler universal handles and are available in 30mm, 45mm, and 60mm staple length options, and in three staple sizes, 2.5mm, 3.5mm, and 4.8mm, to accommodate various tissue thicknesses. Each endoscopic linear cutter stapler reloading unit is single use.	30mm, 45mm, and 60mm staple lengths. 2.5mm, 3.5mm, and 4.8mm staple sizes.	Transferred to ME205
TD003	Purple Surgical Ultimate reloadable linear stapler with indicator - reloads	These are the reloads for the Purple Surgical Ultimate reloadable linear stapler with indicator (Prostheses List application ref. #U5HR-3R6S). The reloads are available in 30mm, 45mm, 60mm, and 90mm staple line lengths for use in different applications. Two staple sizes (3.5mm and 4.8mm) are available to accommodate different tissue thicknesses.	30mm, 45mm, 60mm, and 90mm staple line lengths. 3.5mm and 4.8mm staple sizes.	Transferred to ME204
TD004	Purple Surgical Ultimate haemorrhoidal circular stapler	The Purple Surgical Ultimate haemorrhoidal circular stapler and accessories have application throughout the anal canal to perform the surgical treatment of haemorrhoidal disease. The Purple Surgical Ultimate haemorrhoidal circular stapler places a circular staggered double row of titanium staples in the tissue and resects excess tissue to create a circular anastomosis. Each Purple Surgical Ultimate haemorrhoidal circular stapler comes with 5 accessories: clear circular anal dilator, clear butterfly anal dilator, anal dilator obturator, purse-string suture anoscope, and a suture threader. There are 2 stapler diameters available; 32/33mm and 33/34mm.	32/33mm and 33/34mm.	Transferred to ME200
TD005	Purple Surgical Ultimate reloadable linear stapler with indicator	The Purple Surgical Ultimate linear stapler with indicator can be applied in abdominal, thoracic, and further surgical procedures for transection or resection of tissue. The Purple Surgical Ultimate linear stapler with indicator places a double staggered row of titanium staples and is available in 30mm, 45mm, 60mm, and 90mm staple line lengths for use in different applications. Two staple sizes (3.5mm and 4.8mm) are available to accommodate different tissue thicknesses. Each linear stapler with indicator can be reloaded up to 7 times for a total maximum 8 firings on one patient.	30mm, 45mm, 60mm, and 90mm staple line lengths. 3.5mm and 4.8mm staple sizes .	Transferred to ME203

Billing Code	Product Name	Further Description	Size	Explanation of Deletion
TD006	Purple Surgical Ultimate endoscopic linear cutter stapler reloading units, straight	The Purple Surgical Ultimate endoscopic linear cutter stapler reloading units, straight, are for single patient use only, and are used in conjunction with the Purple Surgical Ultimate endoscopic linear cutter stapler universal handles. A single use reloading unit can be replaced up to 24 times on the same universal handle, during the same procedure. The endoscopic linear cutter stapler straight reloading unit is compatible with the endoscopic linear cutter stapler universal handles and are available in 30mm, 45mm, and 60mm staple length options, and in three staple sizes, 2.5mm, 3.5mm, and 4.8mm, to accommodate various tissue thicknesses. Each endoscopic linear cutter stapler reloading unit is single use.	30mm, 45mm, and 60mm staple lengths. 2.5mm, 3.5mm, and 4.8mm staple sizes.	Transferred to ME201

**Grand Total: 14 Records**